



# Product Catalogue 2024

Bringing the right products and care your way



European Contact Centre & Customer Service Awards

2023

SILVER WINNER

# Laryngectomy Contents

<b>Introduction</b>	<b>4</b>	<b>Accessories</b>	<b>64</b>
What is an HME?	20	Tubeholders and Clips	65
Benefits of using an HME?	21	Protectors	68
About Provox Life	22	Swabs and Brushes	73
		Wipes and Creams	76
		ShowerAids	78
<b>HMEs and Adhesives</b>	<b>18</b>	Foam Covers	80
Provox Life HMEs	24	Tape and Lubricating Jelly	82
Provox Life Adhesives	32	Irrigation Fluid	83
Provox Life Tubes	36	Dressing Packs	84
Provox Life Buttons	40		
Provox HMEs	42		
Provox Adhesives	45		
Provox Tubes	48		
Provox Buttons	51		
Blomsinger System	52		
Spirotect System	60		

# Tracheostomy Contents

<b>Introduction</b>	<b>86</b>	<b>Reversals</b>	<b>118</b>
What is a tracheostomy?	90	Dressing and swabs	119
		Adhesive remover	120
		Skin barrier	120
		Cleaning towels	120
<b>HMEs and Speaking Valves</b>	<b>94</b>		
<b>Accessories</b>	<b>100</b>	<b>Kids</b>	<b>122</b>
Neckbands	101	HMEFs	123
Dressings	103	Neckbands	125
Swabs	106		
Wipes and Creams	108		
Shower Protectors and Tape	110	<b>Jaw Rehabilitation</b>	<b>128</b>
Lubricating Jelly and Irrigation Fluid	111		
Dressing Packs	112		
Protectors	113	<b>Environmental policy</b>	<b>130</b>
		<b>Contacts</b>	<b>131</b>





# Welcome to Atos Care

At Atos Care, our mission is to help make life easier for people living with a neck stoma. We understand that having a laryngectomy or tracheostomy can be life changing, and we are here to support you every step of the way.

In this catalogue you will find information on the services available to you and many of the products we supply. If you can't find what you're looking for, please contact our helpful Customer Care team, who can help you find the products you need if they are available on prescription. We supply products made in-house by Atos (the Provox Life, Provox and Tracoe ranges) and products from a variety of other manufacturers, to give you the best choice and selection.

We are here to give you a helping hand, whatever you may need. If you have any questions or feedback, please do not hesitate to get in touch.

# The Atos Circle of Care

We understand that after your operation, you need more than just products to help you live your life fully. Trust, time and care are more important than ever. The Atos Circle of Care is our comprehensive support service to help you.

**Best Start:** Get off to the best start in life after a laryngectomy or tracheostomy.

The Best Start programme helps you transition from hospital to home, including adjusting to your new condition, understanding your new products and getting in to good routines.

- Enhanced support for your first six months from our Welcome Team of CQC Registered Nurses
- A welcome pack and a welcome call
- Practical tools and equipment to make life easier, including a complimentary care bag containing a range of helpful items
- Regular liaison with clinicians for joined up care



If you would like to find out more about any of our services, please contact our friendly Customer Care team on **0800 783 1659** or email [info@atos-care.co.uk](mailto:info@atos-care.co.uk)



**Connection Hub:** Helping you stay connected right from the beginning of your journey.

Get comfortable with your new life and stay connected to those who know and understand.

- A dedicated Customer Care Representative
- Personalised service - you choose how and when you hear from us
- CQC Registered Nurses to support you in your daily routines, in close partnership and communication with clinicians
- Educational events in the community for people with similar experiences
- Atos MyLife App to provide inspiration and information on living well with a neck stoma



**Care Delivered:** Bringing the right products and care your way.

The right products for your needs, right when you need them.

- Easy and convenient ordering
- Optional monthly reminders
- Electronic Prescription Service
- Rapid, reliable delivery
- Discreet packaging
- Convenience orders by subscription





# How to order

At Atos Care, we are committed to making life easier for people living with a neck stoma. We offer a selection of options for how you can order your supplies; ranging from a monthly service call to a fully automated monthly order.

## 1

**A friendly monthly service call, email or text**

Our helpful Customer Care team can contact you each month to help you place your order and answer any questions you may have.

## 2

**Contact us when you wish to order**

If you prefer to contact us when you need to order, you can contact us using our convenient freephone phone line, or by text, email or through our website, if you prefer.

## 3

**Automated convenience order**

We can set up a convenient repeat order for you, so that your supplies are automatically ordered and sent to you each month.



## Order process explained

- 1 Place an order with our helpful Customer Care team.
- 2 If your order contains new items, we may check suitability with your healthcare professional.
- 3 We contact your GP to request the prescription for your items.\*
- 4 On receiving the prescription\*\*, we send out your supplies directly to your home.

Please allow up to 14 working days from placing your order to receiving your supplies.

Catalogue doesn't contain everything we can provide so if there is anything you need but cannot see then please get in touch with your Customer Care Representative.

\*This applies if your GP agrees to receive prescriptions from a third part (i.e. us). If this is not the case, you may be required to place the request directly with your GP surgery.

\*\*Our Electronic Prescriptions Service (EPS) means your GP can send the prescription directly through us through a secure system, so that you don't have to collect it and send it on to us.



## Dedicated Customer Care team

Our helpful and friendly Customer Care team are here to support you and can help you with understanding, ordering and getting the best out of your products.

You will have a dedicated Customer Care Representative, who will work with you each month:




Helping you to choose the appropriate products for your needs.



Liaising with your clinician to ensure you receive products suitable for your needs.



Establishing a personal relationship to aid you in your journey after surgery.



**“Atos Care are excellent in customer care and the needs of the customer, always professional and a credit to the organisation. Thank you for providing a lovely service and always made me feel welcome.”**

**- Atos Care Service User**

European Contact  
Centre & Customer  
Service Awards

**2023**

**SILVER WINNER**



**93%**

rated the Customer Care team  
overall as good or very good.



## Meet our Atos Care Nurses

To help support you in the community, we provide a complimentary CQC registered Nursing service in the comfort of your home. The service we offer is confidential, discreet and supported by the NHS.\*



\*Please note: At present, our Nursing service is not available in all areas due to local regulations.

Our Atos Care Nursing team is made up of specialist Head and Neck Nurses who can help with a range of important aspects such as:

1

Getting started with Heat and Moisture Exchangers.

2


Guidance on breathing exercises to improve your lung health.

3

Learning to look after your stoma and getting into a good routine.

4

Support you with using your products appropriately.



**Our Atos Care Nursing team work closely with you and your healthcare professionals to understand your individual needs and can provide support via video or phone calls, text messages, emails or a home visit; allowing you choice and flexibility.**

# Events



**MyLife**  
Events

We run a range of events, both online and in-person. Our MyLife events are one of the ways we can help support you to breathe better, whatever you do.

We hold a range of helpful webinars, practical workshops and interactive experiences. Our MyLife events include helpful presentations about relevant topics, as well as a chance to ask questions and hear from other people who know and understand about laryngectomy and tracheostomy.



## Introducing the new MyLife Events



### Webinars

MyLife Webinars are hosted by our expert Nurses and are packed full of information on a specific topic to help improve your quality of life.



### Workshops

MyLife Workshops are hands-on demonstrations held by our specialist team with tips to help you in your everyday life.



### Experiences

MyLife Experiences are online, interactive events for you to meet others and get involved in more activities whether that be a yoga session or a cooking class.



### Community Events

Community events are aimed at bringing communities together and an opportunity for you to meet in-person with those living with a neck stoma in your local area.



### Coffee & Conversations

Come join this virtual event and have a friendly chat over a cuppa. There is no agenda just a chance for you to meet others virtually from the community and ask questions to the Atos team.


Scan to find out more



# Atos MyLife app

Did you know that we have a helpful App, full of informative videos, tutorials and tips about living with a laryngectomy?

Featuring real patient stories as well as presentations from our in-house Speech Therapists and guest experts.



The App is divided into three main sections making it easy for you to find what you are looking for:

## Inspiration

Articles and original videos with patients and clinicians sharing advice and unique stories. Topics include breathing, speaking, care and daily life.

## Exercises

A tool to help you practice voice and breathing exercises through a series of videos offering useful tips to help you overcome any challenges.

## Products

A presentation of the Provox Life products with information to help understand how to use them best within your daily life.





## Award winning app

We are thrilled to announce that our MyLife app won four awards at the 12th annual Lovie awards.

- Gold Lovie Winner: Utilities and Services
- Silver Lovie Winner: Wellness & Lifestyle
- Bronze Lovie Winner: Best practices
- People's Lovie Winner: Best practices

Download the **MyLife app** now



LARYNGECTOMY

# Laryngectomy



# Laryngectomy Contents

<b>Introduction</b>	<b>4</b>	<b>Accessories</b>	<b>64</b>
What is an HME?	20	Tubeholders and Clips	65
Benefits of using an HME?	21	Protectors	68
About Provox Life	22	Swabs and Brushes	73
		Wipes and Creams	76
		ShowerAids	78
<b>HMEs and Adhesives</b>	<b>18</b>	Foam Covers	80
Provox Life HMEs	24	Tape and Lubricating Jelly	82
Provox Life Adhesives	32	Irrigation Fluid	83
Provox Life Tubes	36	Dressing Packs	84
Provox Life Buttons	40		
Provox HMEs	42		
Provox Adhesives	45		
Provox Tubes	48		
Provox Buttons	51		
Blomsinger System	52		
Spirotect System	60		

# What are HMEs?

An HME or Heat and Moisture Exchanger fits onto the end of your laryngectomy tube, laryngectomy button or an adhesive baseplate. You should never put it directly into your stoma. An HME traps heat and moisture from the air you breathe out, which is given back to the air you breathe in.



## Humidification

HMEs humidify (warms up) the air you breathe, which is very important to keep your lungs working properly to help clear mucus away.



## Filtration

HMEs help to filter the air you breathe in by preventing dust, pollen and other particles from entering your lungs.



## Breathing resistance

HMEs provides some breathing resistance which helps you exhale carbon dioxide more efficiently. This is important for your lung health.



# Benefits of using an HME?

Not using HMEs can lead to increased mucus production and coughing. This often can result in increased risk of chest infections, sleep problems, fatigue, less social contact, and an overall reduced quality of life.



**Improve lung health**

Clinical research shows that heat and moisture exchangers (HMEs) have a positive impact on lung health and quality of life after a total laryngectomy, particularly when adhering to HME use 24 hours a day. Using HMEs can also be a more hygienic and easier way of occluding your stoma when speaking with a voice prosthesis.



**Reduce excessive mucus production**



**Reduce coughing**



### Over 1,500 patients in the UK have benefitted from Provox Life

Provox Life™ was developed on the basis of thorough research, and is designed to address the challenges that you experience with HMEs and attachments.

Developed with you in mind, Provox Life™ offers a complete, flexible solution that allows you to experience a personalised routine, no matter where you are on your journey.

The next generation solution consists of:

- Six high performing HMEs that offer improved humidification and breathability.
- Four improved adhesives for a wider range of skin types and stoma contours including a new SecureFit™ coupling.
- Colour coded products and packaging for improved product identification.

**PROVOX**®  
**Life** Breathe better,  
whatever you do



# LARYNGECTOMY

Provox Life™ HMEs are compatible across the whole range of Provox Life™ Adhesives



Provox Life™  
Home HME



Provox Life™  
Go HME



Provox Life™  
Night HME



Provox Life™  
Energy HME



Provox Life™  
Protect HME



Provox Life™  
FreeHands HME



Provox Life™ Stability  
Adhesive



Provox Life™ Night  
Adhesive



Provox Life™ Sensitive  
Adhesive



Provox Life™ Standard  
Adhesive

LARYNGECTOMY

**PROVOX**<sup>®</sup>  
**Life** Breathe better,  
whatever you do

**Provox**<sup>®</sup> **Life**<sup>™</sup>





# LARYNGECTOMY



## Provox® Life™ Home

This HME is designed to improve your lung health when you are taking it easy at home. It features optimised humidification and breathability for healthier lungs.



- Excellent humidification
- Comfortable breathability
- Low profile and discreet look
- Reduces mucus production and coughing
- Easy to press for speaking with a voice prosthesis



## Provox® Life™ Go

This HME provides lower breathing resistance, helping make it easier to breathe when you are out and about and participating in daily activities.



- Good breathability
- Good humidification
- Low-profile and discreet look
- Reduces mucus production and coughing
- Easy to press for speaking with a voice prosthesis

Item No.	Supplier name	Description	Pack size
8311	Atos Medical AB	Provox Life Home HME	30 pc
8310	Atos Medical AB	Provox Life Go HME	30 pc

## LARYNGECTOMY



### Provox® Life™ Night



This HME helps reduce coughing for a more restful nights sleep. It offers our highest humidification which helps reduce mucus production and coughing, allowing you to rest whilst improving lung health.

- Excellent humidification
- Comfortable breathability for nighttime breathing
- Reduces mucus production and coughing
- Soft and smooth for nighttime comfort
- Side openings to prevent occlusion while sleeping



### Provox® Life™ Energy



This HME has good humidification and low breathing resistance, which is perfect for when you are physically active.

- High breathability
- Good humidification
- Sporty look
- Reduces mucus production and coughing
- Easy to press for speaking with a voice prosthesis
- Front vents to prevent accidental occlusion during deep breathing

Item No.	Supplier name	Description	Pack size
8262	Atos Medical AB	Provox Life Night HME	30 pc
8312	Atos Medical AB	Provox Life Energy HME	30 pc

## LARYNGECTOMY



### Provox® Life™ Protect



This HME has a two-way electrostatic filtration that effectively filters bacteria, viruses, dust, and pollen. It also has a comfortable breathing resistance that warms and humidifies the air you breathe.

- A highly effective bacterial and viral filter (greater than or equal to 98%)\*
- Protects you from viruses, bacteria, dust and pollen
- Reduces mucus production and coughing
- Easy to press for speaking with a voice prosthesis

Item No.	Supplier name	Description	Pack size
8313	Atos Medical AB	Provox Life Protect HME	15 pc

**Ready for an upgrade? Text 'LIFE' to 07537 417 928 or call 0800 783 1659**

\*Please note: Since pathogens can enter and leave the human body in other ways (such as the mouth, nose, and eyes), Provox Life™ Protect HME can never guarantee complete protection. Please read the instructions for use for guidance.



## Provox® Life™ FreeHands




Provox® FreeHands HME offers both the freedom to speak hands-free and with finger occlusion. For use with a voice prosthesis. To be used with Provox® FreeHands FlexiVoice™.

- Reduces mucus production and coughing
- Easy removal of HME without the need for a removal aid
- Why not try with the Provox Life™ Stability Adhesive for extra confidence when speaking with Provox Life™ FreeHands

Item No.	Supplier name	Description	Pack size
7440	Atos Medical AB	Provox Life FreeHands HME	30 pc

**PREFERRED PROVOX LIFE™ TO THEIR  
PREVIOUS HMES**

after 12 weeks



Data is on file.

**“I am very happy with the Provox Life™ Stability adhesive and don’t feel nervous it’s going to blow off when voicing.”**

**- Provox Life™ User**

## LARYNGECTOMY



### Provox® FreeHands FlexiVoice™

Gives the ability to speak hands-free. For use with a voice prosthesis.

- Speak hands-free for more natural communication
- Can be covered with a finger if needed
- Quick and easy to set up
- Helpful arch accessory helps keep clothing away
- Available in four strengths to suit your needs
- Compatible with Provox Life™ Adhesives using Provox Life™ FreeHands HME, or with Provox® Adhesives using Provox® FreeHands HMEs

Item No.	Supplier name	Description	Pack size
8161	Atos Medical AB	Provox FreeHands FlexiVoice Light	1 pc
8162	Atos Medical AB	Provox FreeHands FlexiVoice Medium	1 pc
8163	Atos Medical AB	Provox FreeHands FlexiVoice Strong	1 pc
8166	Atos Medical AB	Provox FreeHands FlexiVoice XtraStrong	1 pc

**Your clinician may wish to do a pressure test with you to determine the best FlexiVoice™ for you to use.**

## LARYNGECTOMY



### Provox® FreeHands Support™

To help reduce stoma movement, improve voice quality and prevent loss of baseplate seal while speaking hands-free.

- Helps keep a good seal with your baseplate while speaking hands-free
- Adhesive panel fixes to the chest to reduce baseplate movement whilst speaking
- Three models available to suit different stoma depths

Item No.	Supplier name	Description	Pack size
8021	Atos Medical AB	Provox FreeHands Support Flat	1 pc
8022	Atos Medical AB	Provox FreeHands Support Medium	1 pc
8023	Atos Medical AB	Provox FreeHands Support Deep	1 pc
8024	Atos Medical AB	Provox FreeHands Support Adhesive	15 pcs

## LARYNGECTOMY



### Provox® Life™ Standard Adhesive

An everyday adhesive that provides a strong seal around your stoma. A soft, flexible and low-profile adhesive that forms a strong seal and provides a secure fit.

- Strong seal
- Innovative clover-shaped adhesive for easy application and better fit
- SecureFit™ technology - comfortable low-profile connector
- Compatible with all Provox Life™ HMEs
- Ideal for medium to flat stomas



Item No.	Supplier name	Description	Pack size
7460	Atos Medical AB	Provox Life Standard Adhesive - Round	30 pc
7461	Atos Medical AB	Provox Life Standard Adhesive - Oval	30 pc
7462	Atos Medical AB	Provox Life Standard Adhesive - Plus	30 pc





# LARYNGECTOMY



## Provox® Life™ Sensitive Adhesive

A gentle everyday adhesive for sensitive skin. A comfortable, gentle and low-profile adhesive that forms a good seal for a secure fit.

- Good seal
- Hypoallergenic material for sensitive skin
- Innovative clover-shaped adhesive for easy application and better fit
- SecureFit™ technology - comfortable low-profile connector
- Ideal for sensitive skin
- Compatible with All Provox Life™ HMEs

Item No.	Supplier name	Description	Pack size
7463	Atos Medical AB	Provox Life Sensitive Adhesive - Round	30 pc
7464	Atos Medical AB	Provox Life Sensitive Adhesive - Oval	30 pc
7466	Atos Medical AB	Provox Life Sensitive Adhesive - Plus	30 pc





## Provox® Life™ Stability Adhesive

A firm everyday adhesive that provides a secure attachment, especially if your stoma is deep. It is ideal for those who speak frequently or speak hands-free.



- Strong seal
- Innovative clover-shaped adhesive for easy application and better fit
- SecureFit™ technology - comfortable low-profile connector
- Compatible with all Provox Life™ HMEs
- Ideal for deep stomas or speaking handsfree
- Reinforced centre panel for extra stability

Item No.	Supplier name	Description	Pack size
8263	Atos Medical AB	Provox Life Stability Adhesive	15 pc



## Provox® Life™ Night Adhesive

This adhesive allows your skin to recover while you sleep or during the day. A skin-friendly hydrogel adhesive that soothes and cools the skin.



- Soft and smooth for nighttime comfort
- Hydrogel material helps soothe and cool the skin while you sleep
- Innovative clover-shaped adhesive for easy application and better fit
- SecureFit™ technology - comfortable low-profile connector
- Compatible with all Provox Life™ HMEs
- Ideal for sensitive, irritated skin
- Can be worn during the day when needed

Item No.	Supplier name	Description	Pack size
8261	Atos Medical AB	Provox Life Night Adhesive	15 pc

## LARYNGECTOMY



### Provox® Life™ LaryTube Standard

A flexible silicone tube to attach a Provox Life™ HME or accessory, while making sure the stoma maintains its size and shape.

Can be used with or without a voice prosthesis. To be secured with Provox® TubeHolder or Provox® LaryClips.

Item No.	Supplier name	Description	ID	OD	L	Pack size
7409	Atos Medical AB	Provox Life LaryTube 8/27 - Standard	9.5	12	27	1 pc
7410	Atos Medical AB	Provox Life LaryTube 8/36 - Standard	9.5	12	36	1 pc
7411	Atos Medical AB	Provox Life LaryTube 8/55 - Standard	9.5	12	55	1 pc
7412	Atos Medical AB	Provox Life LaryTube 9/27 - Standard	10.5	13.5	27	1 pc
7413	Atos Medical AB	Provox Life LaryTube 9/36 - Standard	10.5	13.5	36	1 pc
7414	Atos Medical AB	Provox Life LaryTube 9/55 - Standard	10.5	13.5	55	1 pc
7415	Atos Medical AB	Provox Life LaryTube 10/27 - Standard	12	15	27	1 pc
7416	Atos Medical AB	Provox Life LaryTube 10/36 - Standard	12	15	36	1 pc
7417	Atos Medical AB	Provox Life LaryTube 10/55 - Standard	12	15	55	1 pc
7418	Atos Medical AB	Provox Life LaryTube 12/27 - Standard	13.5	17	27	1 pc
7419	Atos Medical AB	Provox Life LaryTube 12/36 - Standard	13.5	17	36	1 pc
7420	Atos Medical AB	Provox Life LaryTube 12/55 - Standard	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

**Please check with your healthcare professional for suitability.**

## LARYNGECTOMY



### Provox® Life™ LaryTube Standard with Ring

A flexible silicone tube to attach a Provox Life™ HME or accessory, while making sure the stoma maintains its size and shape.

Can be used with or without a voice prosthesis. To be secured with Provox® Life™ Adhesives.

Item No.	Supplier name	Description	ID	OD	L	Pack size
7421	Atos Medical AB	Provox Life LaryTube 8/36 - Std with Ring	9.5	12	36	1 pc
7422	Atos Medical AB	Provox Life LaryTube 8/55 - Std with Ring	9.5	12	55	1 pc
7423	Atos Medical AB	Provox Life LaryTube 9/36 - Std with Ring	10.5	13.5	36	1 pc
7424	Atos Medical AB	Provox Life LaryTube 9/55 - Std with Ring	10.5	13.5	55	1 pc
7425	Atos Medical AB	Provox Life LaryTube10/36 - Std with Ring	12	15	36	1 pc
7426	Atos Medical AB	Provox Life LaryTube 10/55 - Std with Ring	12	15	55	1 pc
7427	Atos Medical AB	Provox Life LaryTube 12/36 - Std with Ring	13.5	17	36	1 pc
7428	Atos Medical AB	Provox Life LaryTube 12/55 - Std with Ring	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

**Please check with your healthcare professional for suitability.**

## LARYNGECTOMY



### Provox® Life™ LaryTube Fenestrated

A flexible silicone tube with holes to attach a Provox Life™ HME or accessory, while making sure the stoma maintains its size and shape.

Can be used with or without a voice prosthesis. To be secured with Provox® TubeHolder or Provox® LaryClips.

Item No.	Supplier name	Description	ID	OD	L	Pack size
7429	Atos Medical AB	Provox Life LaryTube 8/36 - Fenestrated	9.5	12	36	1 pc
7430	Atos Medical AB	Provox Life LaryTube 8/55 - Fenestrated	9.5	12	55	1 pc
7431	Atos Medical AB	Provox Life LaryTube 9/36 - Fenestrated	10.5	13.5	36	1 pc
7432	Atos Medical AB	Provox Life LaryTube 9/55 - Fenestrated	10.5	13.5	55	1 pc
7433	Atos Medical AB	Provox Life LaryTube 10/36 - Fenestrated	12	15	36	1 pc
7434	Atos Medical AB	Provox Life LaryTube 10/55 - Fenestrated	12	15	55	1 pc
7435	Atos Medical AB	Provox Life LaryTube 12/36 - Fenestrated	13.5	17	36	1 pc
7436	Atos Medical AB	Provox Life LaryTube 12/55 - Fenestrated	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

Please check with your healthcare professional for suitability.

## LARYNGECTOMY



### Provox® Life™ LaryTube Fenestrated with Ring

A flexible silicone tube with holes to attach a Provox Life™ HME or accessory, while making sure the stoma maintains its size and shape.

Can be used with or without a voice prosthesis. To be secured with Provox® Life™ Adhesives.

Item No.	Supplier name	Description	ID	OD	L	Pack size
8048	Atos Medical AB	Provox Life LaryTube 8/36 - Fen with Ring	9.5	12	36	1 pc
8049	Atos Medical AB	Provox Life LaryTube 8/55 - Fen with Ring	9.5	12	55	1 pc
8050	Atos Medical AB	Provox Life LaryTube 9/36 - Fen with Ring	10.5	13.5	36	1 pc
8051	Atos Medical AB	Provox Life LaryTube 9/55 - Fen with Ring	10.5	13.5	55	1 pc
8052	Atos Medical AB	Provox Life LaryTube 10/36 - Fen with Ring	12	15	36	1 pc
8053	Atos Medical AB	Provox Life LaryTube 10/55 - Fen with Ring	12	15	55	1 pc
8054	Atos Medical AB	Provox Life LaryTube 12/36 - Fen with Ring	13.5	17	36	1 pc
8055	Atos Medical AB	Provox Life LaryTube 12/55 - Fen with Ring	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

Please check with your healthcare professional for suitability.

## LARYNGECTOMY



### Provox® Life™ LaryButton

A soft silicone stoma button to attach a Provox Life™ HME or accessory, while keeping the stoma stable in size.

To be secured with Provox® Tube holder or Provox® LaryClips.

Item No.	Supplier name	Description	OD	L	Pack size
8040	Atos Medical AB	Provox Life LaryButton 12/8	12	8	1 pc
8041	Atos Medical AB	Provox Life LaryButton 12/18	12	18	1 pc
8042	Atos Medical AB	Provox Life LaryButton 14/8	14	8	1 pc
8043	Atos Medical AB	Provox Life LaryButton 14/18	14	18	1 pc
8044	Atos Medical AB	Provox Life LaryButton 16/8	16	8	1 pc
8045	Atos Medical AB	Provox Life LaryButton 16/18	16	18	1 pc
8046	Atos Medical AB	Provox Life LaryButton 18/8	18	8	1 pc
8047	Atos Medical AB	Provox Life LaryButton 18/18	18	18	1 pc


ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

Please check with your healthcare professional for suitability.



## OPTIMISED BREATHABILITY WITH PROVOX LIFE™ HMES

48% reported they never had to take out their HME to catch their breath with Provox Life™.



Data from ReD Associates is on file.

**"I have actually been getting decent sleep since the first night."**

**- Provox Life™ User**



### Provox® XtraMoist™ HME

For excellent humidification when relaxing to help reduce coughing and mucus.

Item No.	Supplier name	Description	Pack size
7290	Atos Medical AB	Provox XtraMoist HME	30 pcs



### Provox® XtraFlow™ HME

For humidification and lower breathing resistance when being active.

Item No.	Supplier name	Description	Pack size
7291	Atos Medical AB	Provox XtraFlow HME	30 pcs

## LARYNGECTOMY



### Provox© Micron™ HME

A highly effective virus and bacteria filter (over 99%\*) with good humidification.

Item No.	Supplier name	Description	Pack size
7248	Atos Medical AB	Provox Micron HME	30 pcs



### Provox© Luna™ HME

A soft HME for excellent humidification when sleeping.

Item No.	Supplier name	Description	Pack size
8013	Atos Medical AB	Provox Luna HME	30 pcs

## LARYNGECTOMY



### Provox® FreeHands HME® Moist

Use with Provox FreeHands FlexiVoice for humidification when relaxing.

Item No.	Supplier name	Description	Pack size
8220	Atos Medical AB	Provox FreeHands HME Moist	30 pcs



### Provox® FreeHands HME® Flow

Use with Provox FreeHands FlexiVoice for a lower breathing resistance when being active.

Item No.	Supplier name	Description	Pack size
8221	Atos Medical AB	Provox FreeHands HME Flow	30 pcs

\*Please note: Since pathogens can enter and leave the human body in other ways (such as the mouth, nose, and eyes), Provox Micro-HME can never guarantee complete protection. Please read the instructions for use for guidance.

## LARYNGECTOMY



### Provox<sup>®</sup> FlexiDerm<sup>™</sup> Adhesive

Gives you a good seal around your stoma.

Item No.	Supplier name	Description	Pack size
7253	Atos Medical AB	Provox FlexiDerm Adhesive - Round	20 pcs
7254	Atos Medical AB	Provox FlexiDerm Adhesive - Oval	20 pcs
7331	Atos Medical AB	Provox FlexiDerm Adhesive - Plus	20 pcs



### Provox<sup>®</sup> OptiDerm<sup>™</sup> Adhesive

An adhesive if your skin feels sore or irritated.

Item No.	Supplier name	Description	Pack size
7255	Atos Medical AB	Provox OptiDerm Adhesive - Round	20 pcs
7256	Atos Medical AB	Provox OptiDerm Adhesive - Oval	20 pcs
7332	Atos Medical AB	Provox OptiDerm Adhesive - Plus	20 pcs

## LARYNGECTOMY



### Provox<sup>®</sup> Stabilibase<sup>™</sup> Adhesive

Gives a superior seal and support when speaking hands-free.

Item No.	Supplier name	Description	Pack size
7289	Atos Medical AB	Provox Stabilibase Adhesive	15 pcs



### Provox<sup>®</sup> Stabilibase OptiDerm<sup>™</sup> Adhesive

When you have a deep stoma and sensitive skin.

Item No.	Supplier name	Description	Pack size
7318	Atos Medical AB	Provox Stabilibase OptiDerm Adhesive	15 pcs

## LARYNGECTOMY



### Provox<sup>®</sup> XtraBase<sup>™</sup> Adhesive

For speaking with finger occlusion or hands-free.

Item No.	Supplier name	Description	Pack size
7265	Atos Medical AB	Provox XtraBase Adhesive	20 pcs



### Provox<sup>®</sup> Luna<sup>™</sup> Adhesive

Designed to help you sleep comfortably and soothe your skin.

Item No.	Supplier name	Description	Pack size
8014	Atos Medical AB	Provox Luna Adhesive	15 pcs



### Provox<sup>®</sup> Adhesive Strip

Protects Provox Life<sup>™</sup> Night Adhesive or Provox<sup>®</sup> Luna Adhesive when showering.

Item No.	Supplier name	Description	Pack size
8015	Atos Medical AB	Provox Adhesive Strip	10 pcs

## LARYNGECTOMY



### Provox® LaryTube Standard

A soft silicone tube used to keep the stoma stable in size.

Can be used with or without a voice prosthesis. To be secured with Provox® TubeHolder or Provox® LaryClips.

Item No.	Supplier name	Description	ID	OD	L	Pack size
7601	Atos Medical AB	Provox LaryTube 8/27 - Standard	9.5	12	27	1 pc
7602	Atos Medical AB	Provox LaryTube 8/36 - Standard	9.5	12	36	1 pc
7603	Atos Medical AB	Provox LaryTube 8/55 - Standard	9.5	12	55	1 pc
7605	Atos Medical AB	Provox LaryTube 9/27 - Standard	10.5	13.5	27	1 pc
7606	Atos Medical AB	Provox LaryTube 9/36 - Standard	10.5	13.5	36	1 pc
7607	Atos Medical AB	Provox LaryTube 9/55 - Standard	10.5	13.5	55	1 pc
7609	Atos Medical AB	Provox LaryTube 10/27 - Standard	12	15	27	1 pc
7610	Atos Medical AB	Provox LaryTube 10/36 - Standard	12	15	36	1 pc
7611	Atos Medical AB	Provox LaryTube 10/55 - Standard	12	15	55	1 pc
7613	Atos Medical AB	Provox LaryTube 12/27 - Standard	13.5	17	27	1 pc
7614	Atos Medical AB	Provox LaryTube 12/36 - Standard	13.5	17	36	1 pc
7615	Atos Medical AB	Provox LaryTube 12/55 - Standard	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

**Please check with your healthcare professional for suitability.**



## LARYNGECTOMY



### Provox® LaryTube Standard with Ring

A soft silicone tube used to keep the stoma stable in size.

Can be used with or without a voice prosthesis. Attaches directly into Provox® adhesive baseplates; no need for neckbands.

Item No.	Supplier name	Description	ID	OD	L	Pack size
7624	Atos Medical AB	Provox LaryTube 8/36 - Std with Ring	9.5	12	36	1 pc
7625	Atos Medical AB	Provox LaryTube 8/55 - Std with Ring	9.5	12	55	1 pc
7626	Atos Medical AB	Provox LaryTube 9/36 - Std with Ring	10.5	13.5	36	1 pc
7627	Atos Medical AB	Provox LaryTube 9/55 - Std with Ring	10.5	13.5	55	1 pc
7628	Atos Medical AB	Provox LaryTube10/36 - Std with Ring	12	15	36	1 pc
7629	Atos Medical AB	Provox LaryTube 10/55 - Std with Ring	12	15	55	1 pc
7630	Atos Medical AB	Provox LaryTube 12/36 - Std with Ring	13.5	17	36	1 pc
7631	Atos Medical AB	Provox LaryTube 12/55 - Std with Ring	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

Please check with your healthcare professional for suitability.

## LARYNGECTOMY



### Provox® LaryTube Fenestrated

A soft silicone tube used to keep the stoma stable in size.

Small holes (fenestrations) are designed for use with a voice prosthesis. To be secured with Provox® TubeHolder or Provox® LaryClips.

Item No.	Supplier name	Description	ID	OD	L	Pack size
7637	Atos Medical AB	Provox LaryTube 8/36 - Fenestrated	9.5	12	36	1 pc
7638	Atos Medical AB	Provox LaryTube 8/55 - Fenestrated	9.5	12	55	1 pc
7640	Atos Medical AB	Provox LaryTube 9/36 - Fenestrated	10.5	13.5	36	1 pc
7641	Atos Medical AB	Provox LaryTube 9/55 - Fenestrated	10.5	13.5	55	1 pc
7643	Atos Medical AB	Provox LaryTube 10/36 - Fenestrated	12	15	36	1 pc
7644	Atos Medical AB	Provox LaryTube 10/55 - Fenestrated	12	15	55	1 pc
7646	Atos Medical AB	Provox LaryTube 12/36 - Fenestrated	13.5	17	36	1 pc
7647	Atos Medical AB	Provox LaryTube 12/55 - Fenestrated	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

Please check with your healthcare professional for suitability.

## LARYNGECTOMY



### Provox® LaryButton

A soft silicone stoma button used to keep the stoma stable in size. To be secured with Provox® TubeHolder or Provox® LaryClips.

Item No.	Supplier name	Description	OD	L	Pack size
7671	Atos Medical AB	Provox LaryButton 12/8	12	8	1 pc
7672	Atos Medical AB	Provox LaryButton 14/8	14	8	1 pc
7673	Atos Medical AB	Provox LaryButton 16/8	16	8	1 pc
7674	Atos Medical AB	Provox LaryButton 18/8	18	8	1 pc
7685	Atos Medical AB	Provox LaryButton 12/18	12	18	1 pc
7686	Atos Medical AB	Provox LaryButton 14/18	14	18	1 pc
7687	Atos Medical AB	Provox LaryButton 16/18	16	18	1 pc
7688	Atos Medical AB	Provox LaryButton 18/18	18	18	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

Please check with your healthcare professional for suitability.

## LARYNGECTOMY

Blom-Singer®



### Blom-Singer® ATSVII Cap & 7 Foam Filters

Twist-on plastic cap with foam disks, for use with the ATSVII Handsfree valve.

Item No.	Supplier name	Description	Pack size
BE1010	Severn Healthcare Technologies	Blom-Singer ATSVII Cap & 7 Foam Filters	1 pc



### Blom-Singer® ATSVII Foam Filters

Replacement foam disks for use with the ATSVII Cap.

Item No.	Supplier name	Description	Pack size
BE1030	Severn Healthcare Technologies	Blom-Singer ATSVII Foam Filters	30 pcs



### Blom-Singer® EasyFlow® HME

Can be pressed down for speaking with a voice prosthesis.

Item No.	Supplier name	Description	Pack size
BE 1055 EF	Severn Healthcare Technologies	Blom-Singer EasyFlow	30 pcs

## LARYNGECTOMY



### Blom-Singer® HME ClassicFlow®

Can be pressed down for speaking with a voice prosthesis.

Item No.	Supplier name	Description	Pack size
BE 1055 EZ	Severn Healthcare Technologies	Blom-Singer HME ClassicFlow	30 pcs



### Blom-Singer® HumidiFilter Holder

Reusable holder that houses the foam filter. Designed for daytime use.

Item No.	Supplier name	Description	Pack size
BE1060	Severn Healthcare Technologies	Blom-Singer HumidiFilter Holder	1 pc



### Blom-Singer® HumidiFilter Holder Foam Filters

A foam disk that can be worn in the re-usable HumidFilter Holder.

Item No.	Supplier name	Description	Pack size
BE 1080	Severn Healthcare Technologies	Blom-Singer HumidiFilter Holder Foam Filters	20 pcs

Please check with your healthcare professional for suitability.

## LARYNGECTOMY



### Blom-Singer® SpeakFree® HME

Single-use adjustable hands-free valve.

Item No.	Supplier name	Description	Pack size
BE 1090 EF	Severn Healthcare Technologies	Blom-Singer SpeakFree HME EasyFlow	30 pcs
BE 1090 EZ	Severn Healthcare Technologies	Blom-Singer SpeakFree HME ClassicFlow	30 pcs



### Blom-Singer® Day and Night® ClassicFlow HME

Touch to speak HME.

Item No.	Supplier name	Description	Pack size
BE 1091	Severn Healthcare Technologies	Blom-Singer Day and Night ClassicFlow	30 pcs



### Blom-Singer® Day and Night® EasyFlow HME

Touch to speak HME.

Item No.	Supplier name	Description	Pack size
BE 1092	Severn Healthcare Technologies	Blom-Singer Day and Night EasyFlow	30 pcs

## LARYNGECTOMY



### Blom-Singer® AccuFit® Adhesive

Flexible, translucent adhesive.

Item No.	Supplier name	Description	Pack size
BE 6082	Severn Healthcare Technologies	Blom-Singer AccuFit Adhesive - Round	30 pcs
BE 6083	Severn Healthcare Technologies	Blom-Singer AccuFit Adhesive - Oval	30 pcs
BE 6084	Severn Healthcare Technologies	Blom-Singer AccuFit Adhesive - Oval Extra	30 pcs



### Blom-Singer® HydroFit® Adhesive

Hydrocolloid adhesive suitable for sensitive skin.

Item No.	Supplier name	Description	Pack size
BE 6085	Severn Healthcare Technologies	Blom-Singer HydroFit Hydrocolloid Adhesive - Round	30 pcs
BE 6086	Severn Healthcare Technologies	Blom-Singer HydroFit Hydrocolloid Adhesive - Oval	30 pcs
BE 6087	Severn Healthcare Technologies	Blom-Singer HydroFit Hydrocolloid Adhesive - Oval Extra	30 pcs

## LARYNGECTOMY



### Blom-Singer® Tracheostoma Valve Housing

A reusable silicone housing for use with Tape Adhesive discs. Fits with the ATSVII or Humidfilter Holder.



Item No.	Supplier name	Description	Pack size
BE6038	Severn Healthcare Technologies	Blom-Singer Tracheostoma Valve Housing - Standard	1 pc
BE6039	Severn Healthcare Technologies	Blom-Singer Tracheostoma Valve Housing - Large	1 pc



## LARYNGECTOMY



### Blom-Singer® Tape Adhesive Discs

Double-sided adhesive tape discs for use with the Tracheostoma Valve Housing.

Item No.	Supplier name	Description	Pack size
BE 6041	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Standard, Regular	30 pcs
BE 6042	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Large, Regular	30 pcs
BE 6043	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Standard, Thin	30 pcs
BE 6044	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Large, Thin	30 pcs
BE6035	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Large, Heavy Duty	30 pcs

## LARYNGECTOMY



### Blom-Singer® Stomasoft® Tube

Soft and flexible laryngectomy tube.

Item No.	Supplier name	Description	ID	OD	L	Pack size
BE 6698	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 8 - Standard	9.5	12	36	1 pc
BE 6699	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 8 - Standard	9.5	12	55	1 pc
BE 6700	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 9 - Standard	10.5	13.5	36	1 pc
BE 6701	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 9 - Standard	10.5	13.5	55	1 pc
BE 6702	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 10 - Standard	12	15	36	1 pc
BE 6703	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 10 - Standard	12	15	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

Please check with your healthcare professional for suitability.

## LARYNGECTOMY



### Blom-Singer® Stomasoft® Tube

Soft and flexible laryngectomy tube with fenestration.

Item No.	Supplier name	Description	ID	OD	L	Pack size
BE 6698F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 8 - Fenestrated	9.5	12	36	1 pc
BE 6699F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 8 - Fenestrated	9.5	12	55	1 pc
BE 6700F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 9 - Fenestrated	10.5	13.5	36	1 pc
BE 6702F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 10 - Fenestrated	12	15	36	1 pc
BE 6703F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 10 - Fenestrated	12	15	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)

Please check with your healthcare professional for suitability.

## LARYNGECTOMY

### SpiroTect©



#### **SpiroTect® Touch2Speak HME**

Easily obtain airflow closure and speech by using a light finger touch.

Item No.	Supplier name	Description	Pack size
PTC1004	Sheffmed Limited	SpiroTect Touch2Speak HME	30 pcs



#### **SpiroTect® Push2Speak HME**

Easily obtain airflow closure and speech by gently pushing the top of the cassette.

Item No.	Supplier name	Description	Pack size
PTC1005	Sheffmed Limited	SpiroTect Push2Speak HME	30 pcs

## LARYNGECTOMY

### SpiroTect® Fabric Adhesive



Fabric dressings that secure HME cassettes.

Item No.	Supplier name	Description	Pack size
PTW100	Sheffmed Limited	Spirotect Fabric Adhesive - Oval	20 pcs
PTW101	Sheffmed Limited	Spirotect Fabric Adhesive - Round	20 pcs

### SpiroTect® Sensitive Adhesive



Made from hydrocolloid that secure HME cassettes.

Item No.	Supplier name	Description	Pack size
PTW200	Sheffmed Limited	Spirotect Sensitive Adhesive - Oval	20 pcs
PTW201	Sheffmed Limited	Spirotect Sensitive Adhesive - Round	20 pcs

## LARYNGECTOMY

### SpiroTect® Clear & Strong Adhesive



Secure HME cassettes at the stoma site.

Item No.	Supplier name	Description	Pack size
PTW300	Sheffmed Limited	SpiroTect Clear & Strong Adhesive - Oval	20 pcs
PTW301	Sheffmed Limited	SpiroTect Clear & Strong Adhesive - Round	20 pcs

### SpiroTect® Extra Sticky Adhesive



Developed with extra adhesion to secure HME cassettes.

Item No.	Supplier name	Description	Pack size
PTW501	Sheffmed Limited	SpiroTect Extra Sticky Adhesive - Round	20 pcs

## LARYNGECTOMY



### **SpiroTect® Secure Adhesive**

Supportive and secure positioning around the stoma site.

Item No.	Supplier name	Description	Pack size
PTW700	Sheffmed Limited	SpiroTect Secure Adhesive - Convex	20 pcs

LARYNGECTOMY

# Accessories





## LARYNGECTOMY

### Accessories - TubeHolders and Clips



#### **Provox® TubeHolder**

This can be used to hold a laryngectomy tube or button in place. The discreet clips fit into the connectors of a Provox® or Provox Life™ LaryTube and both Provox® and Provox Life™ LaryButtons.

- Discreet colour
- Side closure with velcro for adjustable length
- High quality material with slight stretch for comfort
- Compatible with Provox Life™

Item No.	Supplier name	Description	Size	Pack size
7668	Atos Medical AB	Provox TubeHolder	Adjustable	1

## LARYNGECTOMY



### Provox<sup>®</sup> LaryClip

Can be used to hold a laryngectomy tube or button in place. The discreet clips fit into the connectors of a Provox<sup>®</sup> and Provox Life<sup>™</sup> LaryTube or LaryButton.

- Discreet colour
- Adhesive backing attaches securely to skin
- Easy fastenings for flexible positioning
- Compatible with Provox Life<sup>™</sup>

Item No.	Supplier name	Description	Size	Pack size
7669	Atos Medical AB	Provox LaryClip	Adjustable	20 pads/8 clips



### Freevent<sup>®</sup> Neckband One-piece with Twill Ties

- One-piece neckband with twill ties, designed for extra security
- Simply cut and adjust to the exact size needed
- Made of soft foam to help minimise skin irritation
- Rounded edges for comfort

Item No.	Supplier name	Description	Size	Pack size
1661	Atos Medical AB	Freevent Neckband - One-piece	Adult	10

## LARYNGECTOMY



### Freevent® Neckband Two-piece with VELCRO® fastening

- High quality adjustable neckband with VELCRO® fastening designed to provide secure tube positioning
- Simple application with VELCRO® fastening, adjustable to different neck sizes
- Made of soft foam to help minimise skin irritation
- Rounded edges for comfort

Item No.	Supplier name	Description	Size	Pack size
1662	Atos Medical AB	Freevent Neckband - Two-piece	Adult	10



### Insight Tube Holder

Item No.	Supplier name	Description	Size	Pack size
TH100	Insight Medical	Insight Tube Holder - Two-piece	Adult	10

## Accessories - Protectors



### Provox<sup>®</sup> Protector

- Thick, absorbent foam stoma protector
- Ideal for those with lots of secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection, especially in winter

Item No.	Supplier name	Description	Colour	Size	Pack size
7385	Atos Medical AB	Provox Protector	White	Small	10 pcs
7386	Atos Medical AB	Provox Protector	White	Large	10 pcs

## LARYNGECTOMY



### Provox® Protector Slim

- Lightweight, fabric stoma protector
- Easy to breathe through
- Ideal for those with minimal secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection, especially in winter

Item No.	Supplier name	Description	Colour	Size	Pack size
7387	Atos Medical AB	Provox Protector Slim	White	Small	10 pcs
7388	Atos Medical AB	Provox Protector Slim	Blue	Small	10 pcs
7389	Atos Medical AB	Provox Protector Slim	White	Large	10 pcs
7390	Atos Medical AB	Provox Protector Slim	Blue	Large	10 pcs

## LARYNGECTOMY



### Provox® Protector Air

- Light, foam stoma protector
- Easy to breathe through
- Wipe-able - Ideal for those with thick or sticky secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection

Item No.	Supplier name	Description	Colour	Size	Pack size
7391	Atos Medical AB	Provox Protector Air	White	Small	10 pcs
7392	Atos Medical AB	Provox Protector Air	Blue	Small	10 pcs
7393	Atos Medical AB	Provox Protector Air	White	Large	10 pcs
7394	Atos Medical AB	Provox Protector Air	Blue	Large	10 pcs

## LARYNGECTOMY



### Sofshield Protector

Item No.	Supplier name	Description	Colour	Size	Pack size
AS3698	Insight Medical	Sofshield Protector	White	Small	10 pcs
AS3777	Insight Medical	Sofshield Protector	Blue	Small	10 pcs
AS3778	Insight Medical	Sofshield Protector	Beige	Small	10 pcs
AS4056	Insight Medical	Sofshield Protector	Baby blue	Medium	10 pcs
AS3669	Insight Medical	Sofshield Protector	White	Large	10 pcs
AS3764	Insight Medical	Sofshield Protector	Blue	Large	10 pcs
AS3765	Insight Medical	Sofshield Protector	Beige	Large	10 pcs
AS4051	Insight Medical	Sofshield Protector	Black	Large	10 pcs

# LARYNGECTOMY



## Buchanan Protector

Item No.	Supplier name	Description	Colour	Size	Pack size
BUP0002	Sheffmed Limited	Buchanan Protector	White	Small	10 pcs
BUP0001	Sheffmed Limited	Buchanan Protector	White	Large	10 pcs
LABLT2001	Sheffmed Limited	Buchanan Protector Lite	White	Small	10 pcs
LABLT2002	Sheffmed Limited	Buchanan Protector Lite	Blue	Small	10 pcs
LABLT2003	Sheffmed Limited	Buchanan Protector Lite	Beige	Small	10 pcs
LABLT1001	Sheffmed Limited	Buchanan Protector Lite	White	Large	10 pcs
LABLT1002	Sheffmed Limited	Buchanan Protector Lite	Blue	Large	10 pcs
LABLT1003	Sheffmed Limited	Buchanan Protector Lite	Beige	Large	10 pcs



# LARYNGECTOMY

## Accessories - Swabs and Brushes



### Provox<sup>®</sup> Swab

- Strong, sturdy, smooth swab head, designed to clean the surface without scratching
- The head of the swab is made from foam to avoid loose fibres being left behind
- The white swab easily shows blood or coloured mucus
- Robust and flexible handle for ease of movement inside the tube

Item No.	Supplier name	Description	Size	Pack size
8250	Atos Medical AB	Provox Swab Small	4-5 mm	50
8251	Atos Medical AB	Provox Swab Medium	6-8 mm	50
8252	Atos Medical AB	Provox Swab Large	9-10 mm	50
8258	Atos Medical AB	Provox Swab XtraLarge	11 mm-Above	50



### Insight Swab

Item No.	Supplier name	Description	Size	Pack size
TCS/950	Insight Medical	Insight Swab Small	Up to 6 mm	10
TCS/900	Insight Medical	Insight Swab Large	Up to 9 mm	10

## LARYNGECTOMY



### Provox<sup>®</sup> TubeBrush

- Easy to grip handle
- Bristles that extend further down the wire of the brush
- Softer bristles and a soft cotton tip to help reduce the risk of damage to the surface of the tube

Item No.	Supplier name	Description	Size	Pack size
7661	Atos Medical AB	Provox TubeBrush	12 mm	6
7660	Atos Medical AB	Provox TubeBrush	8 mm	6



### Provox<sup>®</sup> Brush

- Rounded tip to protect against accidental injury
- Safety wings to prevent unlikely inhalation of the Provox Brush
- Short handle for easier handling
- Use the Provox Brush to insert the Provox Plug

Item No.	Supplier name	Description	Size	Pack size
7204	Atos Medical AB	Provox Brush	4.5-10 mm	6
7225	Atos Medical AB	Provox Brush XL	12-18 mm	6
8404	Atos Medical AB	Provox Brush Long	4-10 mm	6
8425	Atos Medical AB	Provox Brush XL Long	12-15 mm	6

Please check with your healthcare professional for suitability.

## LARYNGECTOMY



### Blom-Singer® Brush

Item No.	Supplier name	Description	Size	Pack size
IN4064	Severn Healthcare Technologies	Blom-Singer Brush	4-8 mm	3
IN4065	Severn Healthcare Technologies	Blom-Singer Brush	10-16 mm	3
IN4066	Severn Healthcare Technologies	Blom-Singer Brush	18 mm-Above	3

### Insight Brush



Item No.	Supplier name	Description	Size	Pack size
PB/50	Insight Medical	Insight Brush	One size	10

Please check with your healthcare professional for suitability.

## LARYNGECTOMY

### Accessories - Wipes and Creams



#### Provox® Adhesive Remover

- Helps to remove adhesive comfortably and safely
- Reduces risk of skin irritation
- Convenient individual sachets
- No sting formula, safe for use around the stoma

Item No.	Supplier name	Description	Pack size
8012	Atos Medical AB	Provox Adhesive Remover	50



#### Provox® Cleaning Towel

- A convenient and hygienic way to clean around the stoma
- Small pack size, ideal for care bags, handbags and pockets
- No sting formula
- Does not leave any residue on the skin

Item No.	Supplier name	Description	Pack size
7244	Atos Medical AB	Provox Cleaning Towel	10

## LARYNGECTOMY



### Provox<sup>®</sup> Skin Barrier

- Leaves a protective film on the skin to help prevent skin irritation
- Use prior to applying adhesive baseplate
- Convenient individual sachets
- No sting formula, safe for use around the stoma

Item No.	Supplier name	Description	Pack size
8011	Atos Medical AB	Provox Skin Barrier	50



### Spirotect<sup>™</sup> Skin Barrier

Item No.	Supplier name	Description	Pack size
PTA500	Sheffmed Limited	Spirotect Skin Barrier	30



### Spirotect<sup>™</sup> Adhesive Remover

Item No.	Supplier name	Description	Pack size
PTA600	Sheffmed Limited	Spirotect Adhesive Remover	30

## LARYNGECTOMY



### Cavilon© Skin Barrier

Item No.	Supplier name	Description	Pack size
3343P	3M Health Care Ltd	Cavilon Skin Barrier Foam Applicator	5
3391G	3M Health Care Ltd	Cavilon Skin Barrier Cream 28g	1
3392G	3M Health Care Ltd	Cavilon Skin Barrier Cream 92g	1
3344E	3M Health Care Ltd	Cavilon Skin Barrier Wipe	30

### Accessories - Shower Aids



### Provox© Life™ Shower

For use with the Provox Life range.

Item No.	Supplier name	Description	Pack size
8308	Atos Medical AB	Provox Life Shower	1



### Provox© ShowerAid

For use with the Provox range.

Item No.	Supplier name	Description	Pack size
7260	Atos Medical AB	Provox ShowerAid	1

## LARYNGECTOMY



### Provox<sup>®</sup> Luna<sup>™</sup> ShowerAid

For use with Provox Luna Adhesive and Adhesive Strips.

Item No.	Supplier name	Description	Pack size
8016	Atos Medical AB	Provox Luna ShowerAid	1



### Cascade Shower Protector

A flexible, waterproof cover which can be worn over the stoma.

Item No.	Supplier name	Description	Pack size
AS3835	Insight Medical	Cascade Shower Protector	2



### SpiroTect<sup>™</sup> Shower Protector

A flexible, waterproof cover which can be worn over the stoma.

Item No.	Supplier name	Description	Pack size
PTA101	Sheffmed Limited	SpiroTect <sup>™</sup> Shower Protector	2

# LARYNGECTOMY

## Accessories - Foams Covers and Glue



### Provox® Tracheofix

Item No.	Supplier name	Description	Size	Pack size
1429H	Atos Medical AB	Provox Tracheofix 1 strip - Beige	55x60	10
1430H	Atos Medical AB	Provox Tracheofix 1 strip - Beige	70x70	10



Item No.	Supplier name	Description	Size	Pack size
1432H	Atos Medical AB	Provox Tracheofix 2 strips - Beige	55x60	10
1433H	Atos Medical AB	Provox Tracheofix 2 strips - Beige	70x70	10

### Laryngofoam



Item No.	Supplier name	Description	Size	Pack size
LALFF0023	Sheffmed Limited	Laryngofoam - Beige	63x65	30
LALFF0003	Sheffmed Limited	Laryngofoam - Beige	50x62	30



# LARYNGECTOMY



## Purifoam

Item No.	Supplier name	Description	Size	Pack size
CR3853	Insight Medical	Purifoam - Beige	65x65	30



## Provox<sup>®</sup> Adhesive Strip

Protects Provox Life™ Night Adhesive or Provox<sup>®</sup> Luna Adhesive when showering.

Item No.	Supplier name	Description	Size	Pack size
8015	Atos Medical AB	Provox Adhesive Strip		10 pcs



## Provox<sup>®</sup> Silicone Glue

Item No.	Supplier name	Description	Size	Pack size
7720	Atos Medical AB	Provox Silicone Glue	30ml	1

# LARYNGECTOMY

## Accessories - Tape, Lubricating Jelly and Irrigation Fluid



### Insil Tape

Item No.	Supplier name	Description	Size	Pack size
INSIL200	Insight Medical	Insil Soft Silicone Perforated Tape	20mmx3m	1
INSIL400	Insight Medical	Insil Soft Silicone Perforated Tape	40mmx1.5m	1



### Optilube™

Item No.	Supplier name	Description	Size	Pack size
1120	Optimum Medical Solutions	Optilube Lubricating Jelly Sachets	5g	150 sachets
1121	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	42g	1
1122	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	82g	1

## LARYNGECTOMY



### Clinipod<sup>®</sup>

Item No.	Supplier name	Description	Size	Pack size
CLINIPOD	Mayors Healthcare	Clinipod Sodium Chloride 0.9% Irrigation Solution	20ml	25



### Irripod<sup>®</sup>

Item No.	Supplier name	Description	Size	Pack size
IRR24N	C D Medical Ltd	Irripod Sodium Chloride 0.9% Irrigation Solution	20ml	25

## LARYNGECTOMY

### Accessories - Dressing Packs



#### Sterile Dressing Pack Spec 10

Item No.	Supplier name	Description	Pack size
DP SPEC 10	Generic Supplier	Sterile Dressing Pack Spec 10	12 pcs



#### Sterile Dressing Pack Spec 35

Item No.	Supplier name	Description	Pack size
DP SPEC 35	Generic Supplier	Sterile Dressing Pack Spec 35	12 pcs

## LARYNGECTOMY



### Non-Woven Gauze Swabs

Item No.	Supplier name	Description	Pack size
SWAB SPEC 28	Generic Supplier	Non-Woven Gauze Swabs	100

TRACHEOSTOMY

# Tracheostomy



# Tracheostomy Contents

<b>Introduction</b>	<b>86</b>	<b>Reversals</b>	<b>118</b>
What is a tracheostomy?	90	Dressing and swabs	119
		Adhesive remover	120
<b>HMEs and Speaking Valves</b>	<b>94</b>	Skin barrier	120
		Cleaning towels	120
<b>Accessories</b>	<b>100</b>	<b>Kids</b>	<b>122</b>
Neckbands	101	HMEFs	123
Dressings	103	Neckbands	125
Swabs	106		
Wipes and Creams	108	<b>Jaw Rehabilitation</b>	<b>128</b>
Shower Protectors and Tape	110		
Lubricating Jelly and Irrigation Fluid	111	<b>Environmental policy</b>	<b>130</b>
Dressing Packs	112		
Protectors	113	<b>Contacts</b>	<b>131</b>

# TRACHEOSTOMY





# Breathing, speaking and living well

Tracoe® is a tracheostomy brand that includes a range of high quality devices designed to meet high standards for comfort, reliability and durability. HMEs (Heat and Moisture Exchangers) help to restore the lost function of the nose, enabling those with a tracheostomy to breathe better and to speak again.



# What is a tracheostomy?

A tracheostomy can be either temporary or permanent and there are several reasons for having the procedure, all involving restricted airways. After a tracheostomy you will have a tracheostoma, with a tracheostomy tube entering your windpipe (trachea). This will change your way of breathing and speaking.



Before a tracheostomy




After a tracheostomy




## A new way of breathing

Going through a tracheostomy means breathing through a stoma instead of your nose and mouth. Without the nose, air is cooler and less moist than it should be, causing your lungs to produce more mucus.

Heat and Moisture Exchangers (HMEs) will help to restore the lost functions of the nose and mouth to help keep your lungs healthy and working well.



The specially-treated foam inside the HME traps moisture and warmth as you breathe out.



The air passes through the foam inside the HME, warming and moistening the air on its way to the lungs.

## TRACHEOSTOMY



## HME range



### TrachPhone®

The TrachPhone HME is designed to improve your lung health. It can be worn when physically active and when relaxing. It's also suitable to be worn during sleep.



### XtraCare™

The XtraCare HME is designed to provide extra protection by combining an HME with an electrostatic filter. It protects against viruses and bacteria, for when you are away from home.



### DualCare™

The DualCare HME combines a speaking valve with an HME, allowing you to speak hands-free. By twisting the lid of the valve you can switch from speaking mode to HME mode and back.

Apart from the HMEs, there is a variety of additional products that you may need. Tracoe® was developed from valuable experiences given from the many people who have settled back into their daily lives.

## TRACHEOSTOMY



### TrachPhone® HME

TrachPhone is a versatile and lightweight HME that provides good humidification for different situations.

- A multifunctional HME offering effective lung rehabilitation
- Speaking valve function, press the top lid to occlude the stoma for voicing
- Integrated oxygen supply port (4 mm)
- Integrated suctioning port that allows for tracheal suction without removing the HME

Item No.	Supplier name	Description	Pack size
7707	Atos Medical AB	TrachPhone HME	30

## TRACHEOSTOMY



### XtraCare™ HME-F

XtraCare HME provides good humidification and effective filtration, for extra protection.

- Combines an HME with an electrostatic filter, for extra protection
- Protects against airborne particles, including viruses and bacteria, with a filtration efficiency of more than 99%\*
- Can be used with Freevent O<sup>2</sup> Adaptor for administration of supplemental oxygen

Item No.	Supplier name	Description	Pack size
7767	Atos Medical AB	XtraCare HME White	30



### XtraCare™ HME-F Blue

Blue edition is developed to make it easy to see when you are using XtraCare HME-F.

Item No.	Supplier name	Description	Pack size
7768	Atos Medical AB	XtraCare HME Blue	30

## TRACHEOSTOMY



### Freevent® HME

Heat and Moisture Exchanger to be used with DualCare Speaking Valve or DigiTop.

Item No.	Supplier name	Description	Pack size
7742	Atos Medical AB	Freevent HME 15 Regular	30



### O<sup>2</sup> Adaptor

- Adaptor for administration of supplemental oxygen
- To be used with XtraCare

Item No.	Supplier name	Description	Pack size
7769	Atos Medical AB	O2 Adaptor	10



## TRACHEOSTOMY



### DualCare™ Speaking Valve

Combining a speaking valve with the benefits of a Heat and Moisture Exchanger (HME), with a simple twist.

- Unique combination of a speaking valve with a Heat and Moisture Exchanger (HME)
- Change between speaking mode and HME mode with a simple twist
- Hands-free speech in the speaking mode
- Improves or maintains lung health in HME mode

Item No.	Supplier name	Description	Pack size
7744	Atos Medical AB	DualCare Speaking Valve	1



### Freevent® DualCare™ Speaking Valve Blue

Blue edition is developed to make it easy to see when you are using Freevent DualCare Speaking Valve.

Item No.	Supplier name	Description	Pack size
7755	Atos Medical AB	DualCare Speaking Valve Blue	1

## TRACHEOSTOMY

### Trachie-Voice Speaking Valve



Item No.	Supplier name	Description	Pack size
TSV/100	Insight Medical	Trachie-Voice Speaking Valve	2

### Passy Muir® Valve



Item No.	Supplier name	Description	Pack size
TRPMV1002	Passy Muir	Passy Muir Valve 2000 (clear)	1
TRPMV1003	Passy Muir	Passy Muir Valve 2001 (purple)	1
TRPMV0002	Passy Muir	Passy Muir Valve 2001 (Aqua)	1



TRACHEOSTOMY

# Accessories



Accessories - Neckbands

**Freevent® Neckband One-piece with Twill Ties**

- One-piece neckband with twill ties, designed for extra security
- Simply cut and adjust to the exact size needed
- Made of soft foam to help minimise skin irritation
- Rounded edges for comfort



Item No.	Supplier name	Description	Size	Pack size
1661	Atos Medical AB	Freevent Neckband - One-piece	Adult	10

**Freevent® Neckband Two-piece with VELCRO® fastening**

- High quality adjustable neckband with VELCRO® fastening designed to provide secure tube positioning
- Simple application with VELCRO® fastening, adjustable to different neck sizes
- Made of soft foam to help minimise skin irritation
- Rounded edges for comfort



Item No.	Supplier name	Description	Size	Pack size
1662	Atos Medical AB	Freevent Neckband - Two-piece	Adult	10

## TRACHEOSTOMY

### Insight Tube Holder



Item No.	Supplier name	Description	Size	Pack size
TH100	Insight Medical	Insight Tube Holder - Two-piece	Adult	10

### Marpac™ Tube Holder



Item No.	Supplier name	Description	Size	Pack size
201	CCMed	Marpac Tube Holder - One-piece	Adult	10
108D	CCMed	Marpac Tube Holder - One-piece	Adult	25

# TRACHEOSTOMY

## Accessories - Dressings



### Tracoe Purofoam

Tracoe Purofoam dressing is made of a highly absorbent open pore foam. The outer layer is covered with a beige coating which prevents leakage of secretions to the outside and reduces bacterial contamination.

- Allows for quick absorption of secretions
- Non-stick material for easy removal and replacement
- Circular opening and zigzag slit enables easy replacement whilst wearing the tracheostomy tube

Item No.	Supplier name	Description	Size	Pack size
Ref 958	Atos Medical AB	Tracoe Purofoam Large	90mm x 100mm	10
Ref 959	Atos Medical AB	Tracoe Purofoam Small	65mm x 65mm	10

## TRACHEOSTOMY



### Trachi-Dress

- Slimline cotton-type dressing, suitable for those with minimal secretions
- Fenestrated with non-adhesive backing

Item No.	Supplier name	Description	Size	Pack size
TRDRE1001	Atos Medical AB	Trachi-Dress Dressing Small	60mm x 80mm	20
TRDRE1002	Atos Medical AB	Trachi-Dress Dressing Large	80mm x 100mm	20



### Metalline Dressing

- Slimline cotton-type dressing, suitable for those with minimal secretions
- Aluminised backing designed to promote wound healing

Item No.	Supplier name	Description	Size	Pack size
23 094	Lohmann & Rauscher	Metalline Dressing	80mm x 90mm	50



### ActivHeal

- Highly absorbent foam dressing, suitable for those with lots of secretions
- Fenestrated with non-adhesive backing

Item No.	Supplier name	Description	Size	Pack size
10009118	Advanced Medical Solutions	ActivHeal Dressing	100mm x 100mm	10



## TRACHEOSTOMY



### Advazorb (T) Plus

- Highly absorbent foam dressing, suitable for those with lots of secretions
- Fenestrated with non-adhesive backing

Item No.	Supplier name	Description	Size	Pack size
CR3900	Insight Medical	Advazorb (T) Plus Dressing	80 mm x 80 mm	20



### Advazorb (T) Lite

- Slimline foam dressing, suitable for those with minimal secretions
- Fenestrated with non-adhesive backing

Item No.	Supplier name	Description	Size	Pack size
CR4212	Insight Medical	Advazorb (T) Lite Dressing	80 mm x 80 mm	20



### Advadraw (T)

- Suitable for low to medium levels of secretions
- Low adherence - Reduced trauma on removal

Item No.	Supplier name	Description	Size	Pack size
CR4416	Insight Medical	Advadraw (T) Dressing	80 mm x 80 mm	20

# TRACHEOSTOMY

## Accessories - Swabs and Brushes



### Provox<sup>®</sup> Swab

- Strong, sturdy, smooth swab head, designed to clean the surface without scratching
- The head of the swab is made from foam to avoid loose fibres being left behind
- The white swab easily shows blood or coloured mucus
- Robust and flexible handle for ease of movement inside the tube

Item No.	Supplier name	Description	Size	Pack size
8250	Atos Medical AB	Provox Swab Small	4-5 mm	50
8251	Atos Medical AB	Provox Swab Medium	6-8 mm	50
8252	Atos Medical AB	Provox Swab Large	9-10 mm	50
8258	Atos Medical AB	Provox Swab XtraLarge	11 mm-Above	50



### Tracoe<sup>®</sup> Swabs

Item No.	Supplier name	Description	Size	Pack size
TRSWB0150	Atos Medical AB	Trachi-Swab cleaning swab	Small	50
TRSWB0250	Atos Medical AB	Trachi-Swab cleaning swab	Large	50
TRSWB3050	Atos Medical AB	Trachi-Swab cleaning swab	Extra Large	50

# TRACHEOSTOMY



## Insight Swab

Item No.	Supplier name	Description	Size	Pack size
TCS/950	Insight Medical	Insight Swab Small	Up to 6 mm	10
TCS/900	Insight Medical	Insight Swab Large	Up to 9 mm	10

It is recommended that tubes are cleaned with swabs, please check with your healthcare professional for suitability.

## Accessories - Wipes and Creams



### Provox® Cleaning Towel

- A convenient and hygienic way to clean around the stoma
- Small pack size, ideal for care bags, handbags and pockets
- No sting formula
- Does not leave any residue on the skin

Item No.	Supplier name	Description	Pack size
7244	Atos Medical AB	Provox Cleaning Towel	10



### Provox® Skin Barrier

- Leaves a protective film on the skin to help prevent skin irritation
- Use prior to applying adhesive baseplate
- Convenient individual sachets
- No sting formula, safe for use around the stoma

Item No.	Supplier name	Description	Pack size
8011	Atos Medical AB	Provox Skin Barrier	50

## TRACHEOSTOMY



### Cavilon© Skin Barrier

Item No.	Supplier name	Description	Pack size
3343P	3M Health Care Ltd	Cavilon Skin Barrier Foam Applicator	5
3391G	3M Health Care Ltd	Cavilon Skin Barrier Cream 28g	1
3392G	3M Health Care Ltd	Cavilon Skin Barrier Cream 92g	1
3344E	3M Health Care Ltd	Cavilon Skin Barrier Wipe	30

## TRACHEOSTOMY

Accessories - Shower Protector, Tape, Lubricating Jelly and Irrigation Fluid



### Cascade Shower Protector

A flexible, waterproof cover which can be worn over the stoma.

Item No.	Supplier name	Description	Pack size
AS3835	Insight Medical	Cascade Shower Protector	2



### Insil Tape

Item No.	Supplier name	Description	Size	Pack size
INSIL200	Insight Medical	Insil Soft Silicone Perforated Tape	20mmx3m	1
INSIL400	Insight Medical	Insil Soft Silicone Perforated Tape	40mmx1.5m	1

## TRACHEOSTOMY



### Optilube™

Item No.	Supplier name	Description	Size	Pack size
1120	Optimum Medical Solutions	Optilube Lubricating Jelly Sachets	5g	150 sachets
1121	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	42g	1
1122	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	82g	1



### Clinipod©

Item No.	Supplier name	Description	Size	Pack size
CLINIPOD	Mayors Healthcare	Clinipod Sodium Chloride 0.9% Irrigation Solution	20ml	25



### Irripod©

Item No.	Supplier name	Description	Size	Pack size
IRR24N	C D Medical Ltd	Irripod Sodium Chloride 0.9% Irrigation Solution	20ml	25

# TRACHEOSTOMY

## Accessories - Dressing Packs



### Sterile Dressing Pack Spec 10

Item No.	Supplier name	Description	Pack size
DP SPEC 10	Generic Supplier	Sterile Dressing Pack Spec 10	12 pcs



### Sterile Dressing Pack Spec 35

Item No.	Supplier name	Description	Pack size
DP SPEC 35	Generic Supplier	Sterile Dressing Pack Spec 35	12 pcs



### Non-Woven Gauze Swabs

Item No.	Supplier name	Description	Pack size
SWAB SPEC 28	Generic Supplier	Non-Woven Gauze Swabs	100



### Non-Woven Sterile Gauze Swabs

Item No.	Supplier name	Description	Pack size
1860A	Generic Supplier	Non-Woven Sterile Gauze Swabs	200



## Accessories - Protectors



### Provox<sup>®</sup> Protector

- Thick, absorbent foam stoma protector
- Ideal for those with lots of secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection, especially in winter

Item No.	Supplier name	Description	Colour	Size	Pack size
7385	Atos Medical AB	Provox Protector	White	Small	10 pcs
7386	Atos Medical AB	Provox Protector	White	Large	10 pcs



## Provox® Protector Slim

- Lightweight, fabric stoma protector
- Easy to breathe through
- Ideal for those with minimal secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection, especially in winter

Item No.	Supplier name	Description	Colour	Size	Pack size
7387	Atos Medical AB	Provox Protector Slim	White	Small	10 pcs
7388	Atos Medical AB	Provox Protector Slim	Blue	Small	10 pcs
7389	Atos Medical AB	Provox Protector Slim	White	Large	10 pcs
7390	Atos Medical AB	Provox Protector Slim	Blue	Large	10 pcs

## TRACHEOSTOMY



### Provox® Protector Air

- Light, foam stoma protector
- Easy to breathe through
- Wipe-able - Ideal for those with thick or sticky secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection

Item No.	Supplier name	Description	Colour	Size	Pack size
7391	Atos Medical AB	Provox Protector Air	White	Small	10 pcs
7392	Atos Medical AB	Provox Protector Air	Blue	Small	10 pcs
7393	Atos Medical AB	Provox Protector Air	White	Large	10 pcs
7394	Atos Medical AB	Provox Protector Air	Blue	Large	10 pcs

## TRACHEOSTOMY



### Sofshield Protector

Item No.	Supplier name	Description	Colour	Size	Pack size
AS3698	Insight Medical	Sofshield Protector	White	Small	10 pcs
AS3777	Insight Medical	Sofshield Protector	Blue	Small	10 pcs
AS3778	Insight Medical	Sofshield Protector	Beige	Small	10 pcs
AS4056	Insight Medical	Sofshield Protector	Baby blue	Medium	10 pcs
AS3669	Insight Medical	Sofshield Protector	White	Large	10 pcs
AS3764	Insight Medical	Sofshield Protector	Blue	Large	10 pcs
AS3765	Insight Medical	Sofshield Protector	Beige	Large	10 pcs
AS4051	Insight Medical	Sofshield Protector	Black	Large	10 pcs

## TRACHEOSTOMY



### Buchanan Protector

Item No.	Supplier name	Description	Colour	Size	Pack size
BUP0002	Sheffmed Limited	Buchanan Protector	White	Small	10 pcs
BUP0001	Sheffmed Limited	Buchanan Protector	White	Large	10 pcs
LABLT2001	Sheffmed Limited	Buchanan Protector Lite	White	Small	10 pcs
LABLT2002	Sheffmed Limited	Buchanan Protector Lite	Blue	Small	10 pcs
LABLT2003	Sheffmed Limited	Buchanan Protector Lite	Beige	Small	10 pcs
LABLT1001	Sheffmed Limited	Buchanan Protector Lite	White	Large	10 pcs
LABLT1002	Sheffmed Limited	Buchanan Protector Lite	Blue	Large	10 pcs
LABLT1003	Sheffmed Limited	Buchanan Protector Lite	Beige	Large	10 pcs

# Reversals

Once decannulated, you may still need products to care for your neck stoma during the healing process



# TRACHEOSTOMY

## Accessories - Dressing Packs and Swabs



### Antimicrobial Foam Dressing

Item No.	Supplier name	Description	Size	Pack size
66800269	Allevyn	Antimicrobial Foam Dressing	75mm x 75mm	10



### Advazorb Border Lite Foam Dressing

Item No.	Supplier name	Description	Size	Pack size
CR4198	Advancis Medical	Advazorb Border Lite	100mm x 100mm	10



### Non-Woven Sterile Gauze Swabs

Item No.	Supplier name	Description	Pack size
1860A	Generic Supplier	Non-Woven Sterile Gauze Swabs	200

## TRACHEOSTOMY

### Accessories - Wipes and Creams



#### Provox® Adhesive Remover

Item No.	Supplier name	Description	Pack size
8012	Atos Medical AB	Provox Adhesive Remover	50



#### Provox® Skin Barrier

Item No.	Supplier name	Description	Pack size
8011	Atos Medical AB	Provox Skin Barrier	50



#### Provox® Cleaning Towel

Item No.	Supplier name	Description	Pack size
7244	Atos Medical AB	Provox Cleaning Towel	10





# Kids



## TRACHEOSTOMY

### XtraCare™ Mini HMEF



XtraCare Mini provides effective filtration along with good humidification for daily protection.



- Effective filtration of bacteria greater than or equal to 99%\*
- Effective filtration of viruses greater than or equal to 98%\*
- Good humidification
- Compact design, tailored for children with a tracheostomy
- Available in three different colors
- If supplemental oxygen is needed, XtraCare Mini HMEF can be combined with the O<sub>2</sub> Adaptor Mini.



Item No.	Supplier name	Description	Pack size
8004	Atos Medical AB	XtraCare Mini White	30
8005	Atos Medical AB	XtraCare Mini Blue	30
8006	Atos Medical AB	XtraCare Mini Pink	30

\*Please note: Since pathogens can enter and leave the human body in other ways (such as the mouth, nose, and eyes), Freevent® XtraCare Mini HME can never guarantee complete protection. Please read the instructions for use for guidance.

**Please check with your healthcare professional for suitability.**

## TRACHEOSTOMY



### O<sup>2</sup> Adaptor Mini

- Adaptor for administration of supplemental oxygen
- To be used with XtraCare Mini HMEF

Item No.	Supplier name	Description	Pack size
8007	Atos Medical AB	O2 Adaptor Mini	10

Please check with your healthcare professional for suitability.

## TRACHEOSTOMY

### Kids - Neckbands



#### Freevent® Neckband One-piece with Twill Ties

Item No.	Supplier name	Description	Size	Pack size
1651	Atos Medical AB	Freevent Neckband - One-piece	Paediatric	10



#### Freevent® Neckband Two-piece with VELCRO® fastening

Item No.	Supplier name	Description	Size	Pack size
1652	Atos Medical AB	Freevent Neckband - Two-piece	Paediatric	10



#### Insight Tube Holder

Item No.	Supplier name	Description	Size	Pack size
TH150	Insight Medical	Insight Tube Holder - Two-piece	Paediatric	10

## TRACHEOSTOMY



### Marpac™ Tube Holder One-piece with Twill Ties

Item No.	Supplier name	Description	Size	Pack size
106	CC Med	Marpac Tube Holder One-Piece	Neo-natal	10
108D	CCMed	Marpac Tube Holder One-Piece	Paediatric	25



### Marpac™ Tube Holder One-piece with Hybrid Fastener

Item No.	Supplier name	Description	Size	Pack size
107	CC Med	Marpac Tube Holder One-Piece	Paediatric	25
203MARPAC	CC Med	Marpac Tube Holder Two-Piece	Neo-natal/Paediatric	10

Please check with your healthcare professional for suitability.

## TRACHEOSTOMY



# Jaw





## JAW REHABILITATION



### TheraBite®

- For rehabilitation of reduced jaw opening (trismus).
- The TheraBite® Jaw Motion Rehabilitation System™ is a clinically proven, anatomically correct and portable system. Designed to restore mobility after trismus and mandibular hypomobility.
- The TheraBite® Jaw Motion Rehabilitation System™ can also help to reduce muscle pain by reducing joint inflammation. It can prevent reduced muscle strength during radiation treatment (chemotherapy) in the head and neck area.
- The TheraBite Jaw Motion Rehabilitation System is available in both adult and pediatric sizes.

#### Therabite® System Contains:

- 1 x Jaw Mobiliser
- 4 x Bite Pads
- 1 x Hand Aid
- 30 x Range Of Motion Scales
- 1 x Carry Bag
- 1 x Instructions For Use
- 1 x Patient Progress Log

Item No.	Supplier name	Description	Quantity
TH001	Atos Medical AB	TheraBite Jaw Motion Reh. System Adult	1
TH002	Atos Medical AB	TheraBite Jaw Motion Reh. System Pediat.	1
PA001	Atos Medical AB	TheraBite Bite Pad, Regular	4
PA002	Atos Medical AB	TheraBite Bite Pad, Edentulous	4
PA003	Atos Medical AB	TheraBite Bite Pad, Pediatric	4

# Environmental policies



## Fully recyclable packaging

To reduce our carbon footprint we use recyclable materials for every bit of our packaging.



## Delivery using electric vehicles

Most of our deliveries are made using electric vehicles which is being expanded over the next 5 years.



## Responsible forestry

Where possible, our paper and cardboard are sourced from FSC (Forest Stewardship Council) sources.



## Biodegradable materials

Our packaging 'peanuts' are made from biodegradable corn starch and can be dissolved in water.

# Contact us

There are several easy ways to get in touch:



## **Call or text us**

Tel: 0800 783 1659 Text: 0753 741 7928



## **Email us**

[info@atos-care.co.uk](mailto:info@atos-care.co.uk)



## **Visit our website**

[www.atos-care.co.uk](http://www.atos-care.co.uk)



## **Write to us**

Atos Medical, Cartwright House, Riverside Business Park, Tottle Road, Nottingham, NG2 1RT

# Atos

Breathing-Speaking-Living

Atos Care is the distribution arm of Atos Medical, the world-leading supplier of ostomy products for post laryngectomy and tracheostomy patients. Dedicated to making life easier for people living with a neck stoma, Atos Medical has decades of experience, and a deep understanding of what life is like after the operation.

So, although we supply the full range of products, we are much more than just a delivery service. Our expertise and focus mean that we can also provide full professional support and aftercare for you as an individual, wherever you are on your journey. We call this the Atos Circle of Care.

Atos Medical (UK), Cartwright House,  
Tottle Road, Nottingham, NG2 1RT  
Tel: 0800 783 1659 Text: 0753 741 7928  
Email: [info@atos-care.co.uk](mailto:info@atos-care.co.uk)  
Web: [www.atos-care.co.uk](http://www.atos-care.co.uk)