

# Laryngectomy Product Catalogue 2026

Bringing the right products  
and care your way



# Contents

How to order	4
Order process explained	5
Events	6-7
Provox Life HMEs	10-16
Provox Life adhesives	18-21
Provox Life LaryTubes	22-25
Provox Life LaryButton	26
Provox HMEs	28-30
Provox adhesives	31-33
Provox LaryTubes	34-36
Provox LaryButton	37
Accessories	38-50



# Welcome to Atos Care

At Atos Care, our mission is to help make life easier for people living with a neck stoma. We understand that having a laryngectomy can be life changing, and we are here to support you every step of the way.

In this catalogue you will find information on the services available to you and many of the products we supply. If you can't find what you're looking for, please contact our helpful Customer Care team, who can help you find the

products you need if they are available on prescription. We supply products made in-house by Atos (the Provox Life and Provox ranges) and products from a variety of other manufacturers, to give you the best choice and selection.

We are here to give you a helping hand, whatever you may need. If you have any questions or feedback, please do not hesitate to get in touch.

## Contact us

Freephone: **0800 783 1659**

Text: **0753 741 7928**

Email: [info@atos-care.co.uk](mailto:info@atos-care.co.uk)



This catalogue contains Atos laryngectomy products. We can supply all laryngectomy products listed on the NHS Drug Tariff, so if you have additional supplies that you require please contact us to check.

# How to order

At Atos Care, we are committed to making life easier for people living with a neck stoma. We offer a selection of options for how you can order your supplies; ranging from a monthly service call to a fully automated monthly order.



## A monthly service call, email or text

Our helpful Customer Care team can contact you each month to help you place your order and answer any questions you may have.



## Contact us when you wish to order

If you prefer to contact us when you need to order, you can contact us using our convenient freephone phone line, or by text, email or through our website, if you prefer.



## Automated convenience order

We can set up a convenient repeat order for you, so that your supplies are automatically ordered and sent to you each month.

Visit [www.atos-care.co.uk](http://www.atos-care.co.uk) for an easy and convenient way to order products.

# Order process explained



Place an order with our helpful Customer Care team.



If your order contains new items, we may check suitability with your healthcare professional.



We contact your GP to request the prescription for your items.\*



On receiving the prescription\*\*, we send out your supplies directly to your home.

\*This applies if your GP agrees to receive prescriptions from a third party (i.e. us). If this is not the case, you may be required to place the request directly with your GP surgery.

\*\*Our Electronic Prescriptions Service (EPS) means your GP can send the prescription directly through us through a secure system, so that you don't have to collect it and send it on to us.

Please allow up to 14 working days from placing your order to receiving your supplies.

# Events

We run a range of events, both online and in-person. Our MyLife events are one of the ways we can help support you to breathe better, whatever you do.

We hold a range of helpful webinars, practical workshops and interactive experiences. Our MyLife events include helpful presentations about relevant topics, as well as a chance to ask questions and hear from other people who know and understand about laryngectomy.



# Access to experts

Supported events can also include healthcare and other professionals who can share valuable information, tips and best practice that might be hard to find on your own. Past events have included presentations and speakers on themes such as:

- ✓ Looking after your lungs
- ✓ Travel
- ✓ Breathing exercises
- ✓ Dealing with the unexpected




Scan to find out more about our online events.

# MyLife App

Did you know that we have a helpful App, full of informative videos, tutorials and tips about living with a laryngectomy?

Featuring real patient stories as well as presentations from our in-house Speech Therapists, Specialist Nurses and guest experts.

The App is divided into three main sections making it easy for you to find what you are looking for:



## Inspiration

Articles and original videos with patients and clinicians sharing advice and unique stories. Topics include breathing, speaking, care and daily life.



## Excercises

A tool to help you practice voice and breathing exercises through a series of videos offering useful tips to help you overcome any challenges.



## Products

A presentation of the Provox Life products with information to help understand how to use them best within your daily life.



Download the MyLife app now.



Provox Life



## Provox Life Home

This HME is designed to improve your lung health when you are taking it easy at home. It features optimised humidification and breathability for healthier lungs.



- ✓ Excellent humidification
- ✓ Comfortable breathability
- ✓ Low profile and discreet look
- ✓ Reduces mucus production and coughing
- ✓ Easy to press for speaking with a voice prosthesis

## Provox Life Go

This HME provides lower breathing resistance, helping make it easier to breathe when you are out and about and participating in daily activities.



- ✓ Good breathability
- ✓ Good humidification
- ✓ Low-profile and discreet look
- ✓ Reduces mucus production and coughing
- ✓ Easy to press for speaking with a voice prosthesis

Item No.	Description	Pack size
8311	Provox Life Home HME	30 pc
8310	Provox Life Go HME	30 pc

## Provox Life Night



This HME offers our highest humidification which helps reduce mucus production and coughing, allowing you to rest whilst improving lung health.

- ✓ Excellent humidification
- ✓ Comfortable breathability for nighttime breathing
- ✓ Reduces mucus production and coughing
- ✓ Soft and smooth for nighttime comfort
- ✓ Side openings to prevent occlusion while sleeping

## Provox Life Energy



This HME has good humidification and low breathing resistance, which is perfect for when you are physically active.

- ✓ High breathability
- ✓ Good humidification
- ✓ Sporty look
- ✓ Reduces mucus production and coughing
- ✓ Easy to press for speaking with a voice prosthesis
- ✓ Front vents to prevent accidental occlusion during deep breathing

## Provox Life Protect



This HME has a two-way electrostatic filtration that effectively filters bacteria, viruses, dust, and pollen. It also has a comfortable breathing resistance that warms and humidifies the air you breathe.

- ✓ A highly effective bacterial and viral filter (greater than or equal to 98%)\*
- ✓ Reduces mucus production and coughing
- ✓ Protects you from viruses, bacteria, dust and pollen
- ✓ Easy to press for speaking with a voice prosthesis

\*Please note: Since pathogens can enter and leave the human body in other ways (such as the mouth, nose, and eyes), Provox Life Protect HME can never guarantee complete protection. Please read the instructions for use for guidance.

Item No.	Description	Pack size
8262	Provox Life Night HME	30 pc
8312	Provox Life Energy HME	30 pc
8313	Provox Life Protect HME	15 pc


## Provox Life FreeHands



Provox FreeHands HME offers both the freedom to speak hands-free and with finger occlusion. For use with a voice prosthesis. To be used with Provox FreeHands FlexiVoice.

- ✓ Reduces mucus production and coughing
- ✓ Easy removal of HME without the need for a removal aid
- ✓ Why not try with the Provox Life Stability Adhesive for extra confidence when speaking with Provox Life FreeHands

Item No.	Description	Pack size
7440	Provox Life FreeHands HME	30 pc



**“The biggest single thing that has made an improvement to my life is the Night HME.”**

Jon Organ  
Atos Care and Provox Life user



## Provox FreeHands FlexiVoice

Gives the ability to speak hands-free. For use with a voice prosthesis.

- ✓ Speak hands-free for more natural communication
- ✓ Can be covered with a finger if needed
- ✓ Quick and easy to set up
- ✓ Helpful arch accessory helps keep clothing away
- ✓ Available in four strengths to suit your needs
- ✓ Compatible with Provox Life Adhesives using Provox Life FreeHands HME, or with Provox Adhesives using Provox FreeHands HMEs

Item No.	Description	Pack size
8161	Provox FreeHands FlexiVoice Light	1 pc
8162	Provox FreeHands FlexiVoice Medium	1 pc
8163	Provox FreeHands FlexiVoice Strong	1 pc
8166	Provox FreeHands FlexiVoice XtraStrong	1 pc



Your clinician may wish to do a pressure test with you to determine the best FlexiVoice for you to use.




## Provox FreeHands Support

To help reduce stoma movement, improve voice quality and prevent loss of baseplate seal while speaking hands-free.

- ✓ Helps keep a good seal with your baseplate while speaking hands-free
- ✓ Three models available to suit different stoma depths
- ✓ Adhesive panel fixes to the chest to reduce baseplate movement whilst speaking

Item No.	Description	Pack size
8021	Provox FreeHands Support Flat	1 pc
8022	Provox FreeHands Support Medium	1 pc
8023	Provox FreeHands Support Deep	1 pc
8024	Provox FreeHands Support Adhesive	15 pcs



## Provox Life Standard Adhesive

An everyday adhesive that provides a strong seal around your stoma. A soft, flexible and low-profile adhesive that forms a strong seal and provides a secure fit.

- ✓ Strong seal
- ✓ Innovative clover-shaped adhesive for easy application and better fit
- ✓ SecureFit technology - comfortable low-profile connector
- ✓ Compatible with all Provox Life HMEs
- ✓ Ideal for medium to flat stomas

Item No.	Description	Pack size
7460	Provox Life Standard Adhesive - Round	30 pc
7461	Provox Life Standard Adhesive - Oval	30 pc
7462	Provox Life Standard Adhesive - Plus	30 pc





## Provox Life Sensitive Adhesive

A gentle everyday adhesive for sensitive skin. A comfortable, gentle and low-profile adhesive that forms a good seal for a secure fit.

- ✓ Good seal
- ✓ Hypoallergenic material for sensitive skin
- ✓ Innovative clover-shaped adhesive for easy application and better fit
- ✓ SecureFit technology - comfortable low-profile connector
- ✓ Ideal for sensitive skin
- ✓ Compatible with All Provox Life HMEs

Item No.	Description	Pack size
7463	Provox Life Sensitive Adhesive - Round	30 pc
7464	Provox Life Sensitive Adhesive - Oval	30 pc
7466	Provox Life Sensitive Adhesive - Plus	30 pc






## Provox Life Stability Adhesive

A firm everyday adhesive that provides a secure attachment, especially if your stoma is deep. It is ideal for those who speak frequently or speak hands-free.

- ✓ Strong seal
- ✓ Innovative clover-shaped adhesive for easy application and better fit
- ✓ SecureFit technology - comfortable low-profile connector
- ✓ Compatible with all Provox Life HMEs
- ✓ Ideal for deep stomas or speaking handsfree
- ✓ Reinforced centre panel for extra stability

Item No.	Description	Pack size
8263	Provox Life Stability Adhesive	15 pc



## Provox Life Night Adhesive

This adhesive allows your skin to recover while you sleep or during the day. A skin-friendly hydrogel adhesive that soothes and cools the skin.

- ✓ Soft and smooth for nighttime comfort
- ✓ Hydrogel material helps soothe and cool the skin while you sleep
- ✓ Innovative clover-shaped adhesive for easy application and better fit
- ✓ SecureFit technology - comfortable low-profile connector
- ✓ Compatible with all Provox Life HMEs
- ✓ Ideal for sensitive, irritated skin
- ✓ Can be worn during the day when needed

Item No.	Description	Pack size
8261	Provox Life Night Adhesive	15 pc



## Provox Life LaryTube Standard

A flexible silicone tube to attach a Provox Life HME or accessory, while making sure the stoma maintains its size and shape.

- ✓ Can be used with or without a voice prosthesis. To be secured with Provox TubeHolder or Provox LaryClips

Item No.	Description	ID	OD	L	Pack size
7409	Provox Life LaryTube 8/27 - Standard	9.5	12	27	1 pc
7410	Provox Life LaryTube 8/36 - Standard	9.5	12	36	1 pc
7411	Provox Life LaryTube 8/55 - Standard	9.5	12	55	1 pc
7412	Provox Life LaryTube 9/27 - Standard	10.5	13.5	27	1 pc
7413	Provox Life LaryTube 9/36 - Standard	10.5	13.5	36	1 pc
7414	Provox Life LaryTube 9/55 - Standard	10.5	13.5	55	1 pc
7415	Provox Life LaryTube 10/27 - Standard	12	15	27	1 pc
7416	Provox Life LaryTube 10/36 - Standard	12	15	36	1 pc
7417	Provox Life LaryTube 10/55 - Standard	12	15	55	1 pc
7418	Provox Life LaryTube 12/27 - Standard	13.5	17	27	1 pc
7419	Provox Life LaryTube 12/36 - Standard	13.5	17	36	1 pc
7420	Provox Life LaryTube 12/55 - Standard	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.



## Provox Life LaryTube Standard with Ring

A flexible silicone tube to attach a Provox Life HME or accessory, while making sure the stoma maintains its size and shape.

- ✓ Can be used with or without a voice prosthesis. To be secured with Provox Life Adhesives

Item No.	Description	ID	OD	L	Pack size
7421	Provox Life LaryTube 8/36 - Std with Ring	9.5	12	36	1 pc
7422	Provox Life LaryTube 8/55 - Std with Ring	9.5	12	55	1 pc
7423	Provox Life LaryTube 9/36 - Std with Ring	10.5	13.5	36	1 pc
7424	Provox Life LaryTube 9/55 - Std with Ring	10.5	13.5	55	1 pc
7425	Provox Life LaryTube10/36 - Std with Ring	12	15	36	1 pc
7426	Provox Life LaryTube 10/55 - Std with Ring	12	15	55	1 pc
7427	Provox Life LaryTube 12/36 - Std with Ring	13.5	17	36	1 pc
7428	Provox Life LaryTube 12/55 - Std with Ring	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.



## Provox Life LaryTube Fenestrated

A flexible silicone tube with holes to attach a Provox Life HME or accessory, while making sure the stoma maintains its size and shape.

- ✓ Can be used with or without a voice prosthesis. To be secured with Provox TubeHolder or Provox LaryClips

Item No.	Description	ID	OD	L	Pack size
7429	Provox Life LaryTube 8/36 - Fenestrated	9.5	12	36	1 pc
7430	Provox Life LaryTube 8/55 - Fenestrated	9.5	12	55	1 pc
7431	Provox Life LaryTube 9/36 - Fenestrated	10.5	13.5	36	1 pc
7432	Provox Life LaryTube 9/55 - Fenestrated	10.5	13.5	55	1 pc
7433	Provox Life LaryTube 10/36 - Fenestrated	12	15	36	1 pc
7434	Provox Life LaryTube 10/55 - Fenestrated	12	15	55	1 pc
7435	Provox Life LaryTube 12/36 - Fenestrated	13.5	17	36	1 pc
7436	Provox Life LaryTube 12/55 - Fenestrated	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.



## Provox Life LaryTube Fenestrated with Ring

A flexible silicone tube with holes to attach a Provox Life HME or accessory, while making sure the stoma maintains its size and shape.

- ✓ Can be used with or without a voice prosthesis. To be secured with Provox Life Adhesives

Item No.	Description	ID	OD	L	Pack size
8048	Provox Life LaryTube 8/36 - Fen with Ring	9.5	12	36	1 pc
8049	Provox Life LaryTube 8/55 - Fen with Ring	9.5	12	55	1 pc
8050	Provox Life LaryTube 9/36 - Fen with Ring	10.5	13.5	36	1 pc
8051	Provox Life LaryTube 9/55 - Fen with Ring	10.5	13.5	55	1 pc
8052	Provox Life LaryTube 10/36 - Fen with Ring	12	15	36	1 pc
8053	Provox Life LaryTube 10/55 - Fen with Ring	12	15	55	1 pc
8054	Provox Life LaryTube 12/36 - Fen with Ring	13.5	17	36	1 pc
8055	Provox Life LaryTube 12/55 - Fen with Ring	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.



## Provox Life LaryButton

A soft silicone stoma button to attach a Provox Life HME or accessory, while keeping the stoma stable in size.


✓ To be secured with Provox Tube holder or Provox LaryClips

Item No.	Description	OD	L	Pack size
8040	Provox Life LaryButton 12/8	12	8	1 pc
8041	Provox Life LaryButton 12/18	12	18	1 pc
8042	Provox Life LaryButton 14/8	14	8	1 pc
8043	Provox Life LaryButton 14/18	14	18	1 pc
8044	Provox Life LaryButton 16/8	16	8	1 pc
8045	Provox Life LaryButton 16/18	16	18	1 pc
8046	Provox Life LaryButton 18/8	18	8	1 pc
8047	Provox Life LaryButton 18/18	18	18	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.

A photograph of an elderly man with a white beard and glasses, wearing a white hospital gown, sitting up in a hospital bed. He is holding a smartphone in his hands and looking at the screen. The room is dimly lit, with a bedside lamp glowing and a potted plant on a table next to the bed.

“Ordering my supplies is always easy as I receive a call every month from a dedicated customer service representative who knows me well and exactly what I require on a regular basis.

The service I receive from the customer service department is always friendly and helpful.”

Provox Life user



## Provox XtraMoist HME

For excellent humidification when relaxing to help reduce coughing and mucus.

Item No.	Description	Pack size
7290	Provox XtraMoist HME	30 pcs



## Provox XtraFlow HME

For humidification and lower breathing resistance when being active.

Item No.	Description	Pack size
7291	Provox XtraFlow HME	30 pcs



## Provox Micron HME

A highly effective virus and bacteria filter (over 99%\*) with good humidification.

Item No.	Description	Pack size
7248	Provox Micron HME	30 pcs



## Provox Luna HME

A soft HME for excellent humidification when sleeping.

✓ This HME should be used with Provox Luna adhesive

Item No.	Description	Pack size
8013	Provox Luna HME	30 pcs



## Provox FreeHands HME Moist

Use with Provox FreeHands FlexiVoice for humidification when relaxing.

Item No.	Description	Pack size
8220	Provox FreeHands HME Moist	30 pcs



## Provox FreeHands HME Flow

Use with Provox FreeHands FlexiVoice for a lower breathing resistance when being active.

Item No.	Description	Pack size
8221	Provox FreeHands HME Flow	30 pcs



## Provox FlexiDerm Adhesive

Gives you a good seal around your stoma.

Item No.	Description	Pack size
7253	Provox FlexiDerm Adhesive - Round	20 pcs
7254	Provox FlexiDerm Adhesive - Oval	20 pcs
7331	Provox FlexiDerm Adhesive - Plus	20 pcs



## Provox OptiDerm Adhesive

An adhesive if your skin feels sore or irritated.

Item No.	Description	Pack size
7255	Provox OptiDerm Adhesive - Round	20 pcs
7256	Provox OptiDerm Adhesive - Oval	20 pcs
7332	Provox OptiDerm Adhesive - Plus	20 pcs



## Provox StabiliBase Adhesive

Gives a superior seal and support when speaking hands-free.

- ✓ Acrylic adhesive with firm concave base, providing support and strong adhesion

Item No.	Description	Pack size
7289	Provox StabiliBase Adhesive	15 pcs



## Provox StabiliBase OptiDerm Adhesive

When you have a deep stoma and sensitive skin.

Item No.	Description	Pack size
7318	Provox StabiliBase OptiDerm Adhesive	15 pcs



## Provox XtraBase Adhesive

For speaking with finger occlusion or hands-free.

- ✓ Strong, acrylic adhesive with reinforced concave base for extra support

Item No.	Description	Pack size
7265	Provox XtraBase Adhesive	20 pcs



## Provox Luna Adhesive

Designed to help you sleep comfortably and soothe your skin.

Item No.	Description	Pack size
8014	Provox Luna Adhesive	15 pcs



## Provox Adhesive Strip

Protects Provox Life Night Adhesive or Provox Luna Adhesive when showering.

Item No.	Description	Pack size
8015	Provox Adhesive Strip	10 pcs



## Provox LaryTube Standard

A soft silicone tube used to keep the stoma stable in size.

- ✓ Can be used with or without a voice prosthesis. To be secured with Provox TubeHolder or Provox LaryClips

Item No.	Description	ID	OD	L	Pack size
7601	Provox LaryTube 8/27 - Standard	9.5	12	27	1 pc
7602	Provox LaryTube 8/36 - Standard	9.5	12	36	1 pc
7603	Provox LaryTube 8/55 - Standard	9.5	12	55	1 pc
7605	Provox LaryTube 9/27 - Standard	10.5	13.5	27	1 pc
7606	Provox LaryTube 9/36 - Standard	10.5	13.5	36	1 pc
7607	Provox LaryTube 9/55 - Standard	10.5	13.5	55	1 pc
7609	Provox LaryTube 10/27 - Standard	12	15	27	1 pc
7610	Provox LaryTube 10/36 - Standard	12	15	36	1 pc
7611	Provox LaryTube 10/55 - Standard	12	15	55	1 pc
7613	Provox LaryTube 12/27 - Standard	13.5	17	27	1 pc
7614	Provox LaryTube 12/36 - Standard	13.5	17	36	1 pc
7615	Provox LaryTube 12/55 - Standard	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.



## Provox LaryTube Standard with Ring

A soft silicone tube used to keep the stoma stable in size.

- ✓ Can be used with or without a voice prosthesis. Attaches directly into Provox adhesive baseplates; no need for neckbands

Item No.	Description	ID	OD	L	Pack size
7624	Provox LaryTube 8/36 - Std with Ring	9.5	12	36	1 pc
7625	Provox LaryTube 8/55 - Std with Ring	9.5	12	55	1 pc
7626	Provox LaryTube 9/36 - Std with Ring	10.5	13.5	36	1 pc
7627	Provox LaryTube 9/55 - Std with Ring	10.5	13.5	55	1 pc
7628	Provox LaryTube10/36 - Std with Ring	12	15	36	1 pc
7629	Provox LaryTube 10/55 - Std with Ring	12	15	55	1 pc
7630	Provox LaryTube 12/36 - Std with Ring	13.5	17	36	1 pc
7631	Provox LaryTube 12/55 - Std with Ring	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.



## Provox LaryTube Fenestrated

A soft silicone tube used to keep the stoma stable in size.

- ✓ Small holes (fenestrations) are designed for use with a voice prosthesis. To be secured with Provox TubeHolder or Provox LaryClips

Item No.	Description	ID	OD	L	Pack size
7637	Provox LaryTube 8/36 - Fenestrated	9.5	12	36	1 pc
7638	Provox LaryTube 8/55 - Fenestrated	9.5	12	55	1 pc
7640	Provox LaryTube 9/36 - Fenestrated	10.5	13.5	36	1 pc
7641	Provox LaryTube 9/55 - Fenestrated	10.5	13.5	55	1 pc
7643	Provox LaryTube 10/36 - Fenestrated	12	15	36	1 pc
7644	Provox LaryTube 10/55 - Fenestrated	12	15	55	1 pc
7646	Provox LaryTube 12/36 - Fenestrated	13.5	17	36	1 pc
7647	Provox LaryTube 12/55 - Fenestrated	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.



## Provox LaryButton

A soft silicone stoma button used to keep the stoma stable in size. To be secured with Provox TubeHolder or Provox LaryClips.

Item No.	Description	OD	L	Pack size
7671	Provox LaryButton 12/8	12	8	1 pc
7672	Provox LaryButton 14/8	14	8	1 pc
7673	Provox LaryButton 16/8	16	8	1 pc
7674	Provox LaryButton 18/8	18	8	1 pc
7685	Provox LaryButton 12/18	12	18	1 pc
7686	Provox LaryButton 14/18	14	18	1 pc
7687	Provox LaryButton 16/18	16	18	1 pc
7688	Provox LaryButton 18/18	18	18	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm)  
L = Length (mm)



Please check with your healthcare professional for suitability.

# Accessories





## Provox TubeHolder

This can be used to hold a laryngectomy tube or button in place. The discreet clips fit into the connectors of a Provox or Provox Life LaryTube and both Provox and Provox Life LaryButtons.

- ✓ Side closure with velcro for adjustable length
- ✓ High quality material with slight stretch for comfort
- ✓ Compatible with Provox Life

Item No.	Description	Size	Pack size
7668	Provox TubeHolder	Adjustable	1



## Provox LaryClip

Can be used to hold a laryngectomy tube or button in place. The discreet clips fit into the connectors of a Provox and Provox Life LaryTube or LaryButton.

- ✓ Adhesive backing attaches securely to skin
- ✓ Easy fastenings for flexible positioning
- ✓ Compatible with Provox Life

Item No.	Description	Size	Pack size
7669	Provox LaryClip	Adjustable	20 pads/8 clips



## Freeevent Neckband One-piece with Twill Ties

- ✓ One-piece neckband with twill ties, designed for extra security
- ✓ Simply cut and adjust to the exact size needed
- ✓ Made of soft foam to help minimise skin irritation
- ✓ Rounded edges for comfort

Item No.	Description	Size	Pack size
1661	Freeevent Neckband - One-piece	Upto 550mm	10



## Freeevent Neckband Two-piece with VELCRO fastening

- ✓ High quality adjustable neckband with VELCRO fastening designed to provide secure tube positioning
- ✓ Simple application with VELCRO fastening, adjustable to different neck sizes
- ✓ Made of soft foam to help minimise skin irritation
- ✓ Rounded edges for comfort

Item No.	Description	Size	Pack size
1662	Freeevent Neckband - Two-piece	Upto 550mm	10



## Provox Protector

- ✓ Thick, absorbent foam stoma protector
- ✓ Ideal for those with lots of secretions
- ✓ Comfortable and discreet
- ✓ Soft-lock mechanism for comfort
- ✓ Hand washable up to 3 times
- ✓ Suitable to be worn over an HME for additional protection, especially in winter



## Provox Protector Slim

- ✓ Lightweight, fabric stoma protector
- ✓ Easy to breathe through
- ✓ Ideal for those with minimal secretions
- ✓ Comfortable and discreet
- ✓ Soft-lock mechanism for comfort
- ✓ Hand washable up to 3 times
- ✓ Suitable to be worn over an HME for additional protection, especially in winter



## Provox Protector Air

- ✓ Light, foam stoma protector
- ✓ Easy to breathe through
- ✓ Wipe-able - Ideal for those with thick or sticky secretions
- ✓ Comfortable and discreet
- ✓ Soft-lock mechanism for comfort
- ✓ Hand washable up to 3 times
- ✓ Suitable to be worn over an HME for additional protection

Item No.	Description	Colour	Size	Pack size
7385	Provox Protector	White	Small	10 pcs
7386	Provox Protector	White	Large	10 pcs
7387	Provox Protector Slim	White	Small	10 pcs
7388	Provox Protector Slim	Blue	Small	10 pcs
7389	Provox Protector Slim	White	Large	10 pcs
7390	Provox Protector Slim	Blue	Large	10 pcs
7391	Provox Protector Air	White	Small	10 pcs
7392	Provox Protector Air	Blue	Small	10 pcs
7393	Provox Protector Air	White	Large	10 pcs
7394	Provox Protector Air	Blue	Large	10 pcs



## Provox Swab

- ✓ Strong, sturdy, smooth swab head, designed to clean the surface without scratching
- ✓ The head of the swab is made from foam to avoid loose fibres being left behind
- ✓ The white swab easily shows blood or coloured mucus
- ✓ Robust and flexible handle for ease of movement inside the tube

Item No.	Description	Size	Pack size
8250	Provox Swab Small	4-5 mm	50
8251	Provox Swab Medium	6-8 mm	50
8252	Provox Swab Large	9-10 mm	50
8258	Provox Swab XtraLarge	11 mm-Above	50



## Tracoe Swabs

Item No.	Description	Size	Pack size
TRSWB0150	Trachi-Swab cleaning swab	3.5-5mm	50
TRSWB0250	Trachi-Swab cleaning swab	5-7.5mm	50
TRSWB3050	Trachi-Swab cleaning swab	8mm+	50



## Provox TubeBrush

- ✓ Easy to grip handle
- ✓ Bristles that extend further down the wire of the brush
- ✓ Softer bristles and a soft cotton tip to help reduce the risk of damage to the surface of the tube

Item No.	Description	Size	Pack size
7661	Provox TubeBrush	12 mm	6
7660	Provox TubeBrush	8 mm	6



## Provox Voice Prosthesis Brush

- ✓ Rounded tip to protect against accidental injury
- ✓ Safety wings to prevent unlikely inhalation of the Provox Brush
- ✓ Short handle for easier handling

Item No.	Description	Size	Pack size
7204	Provox Brush	4.5-10 mm	6
7225	Provox Brush XL	12-18 mm	6
8404	Provox Brush Long	4-10 mm	6
8425	Provox Brush XL Long	12-15 mm	6

## Provox Adhesive Remover



- ✓ Helps to remove adhesive comfortably and safely
- ✓ Reduces risk of skin irritation
- ✓ Convenient individual sachets
- ✓ No sting formula, safe for use around the stoma

Item No.	Description	Pack size
8012	Provox Adhesive Remover	50



## Provox Cleaning Towel

- ✓ A convenient and hygienic way to clean around the stoma
- ✓ Small pack size, ideal for care bags, handbags and pockets
- ✓ No sting formula
- ✓ Does not leave any residue on the skin

Item No.	Description	Pack size
7244	Provox Cleaning Towel	10



## Provox Skin Barrier

- ✓ Leaves a protective film on the skin to help prevent skin irritation
- ✓ Use prior to applying adhesive baseplate
- ✓ Convenient individual sachets
- ✓ No sting formula, safe for use around the stoma

Item No.	Description	Pack size
8011	Provox Skin Barrier	50



## Provox Life Shower

For use with the Provox Life range.

Item No.	Description	Pack size
8308	Provox Life Shower	1



## Provox ShowerAid

For use with the Provox range.

Item No.	Description	Pack size
7260	Provox ShowerAid	1



## Provox Luna ShowerAid

For use with Provox Luna Adhesive and Adhesive Strips.

Item No.	Description	Pack size
8016	Provox Luna ShowerAid	1



## Provox Tracheofix

Item No.	Description	Size	Pack size
1429H	Provox Tracheofix 1 strip - Beige	55x60	10
1430H	Provox Tracheofix 1 strip - Beige	70x70	10



Item No.	Description	Size	Pack size
1432H	Provox Tracheofix 2 strips - Beige	55x60	10
1433H	Provox Tracheofix 2 strips - Beige	70x70	10



## Provox Adhesive Strip

Protects Provox Life™ Night Adhesive or Provox Luna Adhesive when showering.

Item No.	Description	Pack size
8015	Provox Adhesive Strip	10 pcs



## Provox Silicone Glue

Item No.	Description	Size	Pack size
7720	Provox Silicone Glue	30ml	1

# Environmental policies



## Fully recyclable packaging

To reduce our carbon footprint we use recyclable materials for every bit of our packaging.



## Responsible forestry

Where possible, our paper and cardboard are sourced from FSC (Forest Stewardship Council) sources.



## Delivery using electric vehicles

Most of our deliveries are made using electric vehicles which is being expanded over the next 5 years.



## Biodegradable materials

Our packaging 'peanuts' are made from biodegradable corn starch and can be dissolved in water.

Atos Care is the distribution arm of Atos Medical, the world-leading supplier of ostomy products for post laryngectomy and tracheostomy patients. Dedicated to making life easier for people living with a neck stoma, Atos Medical has decades of experience, and a deep understanding of what life is like after the operation.

So, although we supply the full range of products, we are much more than just a delivery service. Our expertise and focus mean that we can also provide full professional support and aftercare for you as an individual, wherever you are on your journey. We call this the Atos Circle of Care.

Atos Medical (UK), Cartwright House,  
Tottle Road, Nottingham, NG2 1RT  
Tel: 0800 783 1659  
Text: 0753 741 7928  
Email: [info@atos-care.co.uk](mailto:info@atos-care.co.uk)  
Web: [www.atos-care.co.uk](http://www.atos-care.co.uk)